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GST Gene Fusion System
Handbook
18-1157-58

Affinity Chromatography
Principles and Methods
18-1022-29

Antibody Purification
Handbook
18-1037-46

Cell Separation Media
Methodology and applications
18-1115-69

Isolation of mononuclear cells
Methodology and applications
18-1152-69

**Ion Exchange Chromatography
& Chromatofocusing**
Principles and Methods
11-0004-21

**Purifying
Challenging Proteins**
Principles and Methods
28-9095-31

Gel Filtration
Principles and Methods
18-1022-18

**Recombinant Protein
Purification Handbook**
Principles and Methods
18-1142-75

Protein Purification
Handbook
18-1132-29

**Hydrophobic Interaction and
Reversed Phase Chromatography**
Principles and Methods
11-0012-69

**2-D Electrophoresis
using immobilized pH gradients**
Principles and Methods
80-6429-60

Microcarrier Cell Culture
Principles and Methods
18-1140-62

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1. Isolation of mononuclear cells and granulocytes from human blood. (Paper IV). Bøyum, A. *Scand J Clin Lab Invest 21 Suppl*, **97**, 77-89 (1968).
2. Isolation of leucocytes from human blood - further observations. (Paper II). Bøyum, A. *Scand J Clin Lab Invest 21 Suppl*, **97**, 31-50 (1968).
3. Isolation of lymphocytes, granulocytes and macrophages. Bøyum, A. *Scand J Immunol 5 Suppl*, **5**, 9-15 (1976).
4. Gel formation with leucocytes and heparin. Almeida, A.P., Beaven, M.A. *Life Sci*, **26**, 549-555 (1980).
5. Tissue typing using a routine one-step lymphocyte separation procedure. Harris, R., Ukaejiofo, E.O. *Brit J Haematol*, **18**, 229-235 (1970).
6. Altered lymphocyte markers and blastogenic responses associated with 24 hour delay in processing of blood samples. Kaplan, J., Nolan, D., Ree, A. *J Immunol Methods*, **50**, 187-191 (1982).
7. Purification of lymphocytes and platelets by gradient centrifugation. Perper, R.J., Zee, T.W., Mickelson, M.M. *J Lab Clin Med*, **72**, 842-848 (1968).
8. Platelet aggregation technique used in the preparation of lymphocyte suspensions. Vives, J., Parra, M., Castillo, R. *Tissue Antigens*, **1**, 276-278 (1971).
9. Quantitation of F_c receptors and surface immunoglobulin is affected by cell isolation procedures using Plasmagel and Ficoll-Hypaque. Alexander, E.L., Titus, J.A., Sega, D.M. *J Immunol Methods*, **22**, 263-272 (1978).
10. Analysis of the lymphocyte distribution during Isopaque-Ficoll isolation of mononuclear cells from human peripheral blood. Hokland, P., Heron, I. *J Immunol Methods*, **32**, 31-39 (1980).
11. The Isopaque-Ficoll method re-evaluated: Selective loss of autologous rosette-forming lymphocytes during isolation of mononuclear cells from human peripheral blood. Hokland, P., Heron, I. *Scand J Immunol*, **11**, 353-356 (1980).
12. Reactivity in mixed cultures of mononuclear leucocytes separated on Ficoll-Hypaque. Bain, B., Pshyk, K. *Proceedings 7th Leucocyte Culture Conference* (Ed. Daguillard, F.) *Academic Press*, New York, 29-37 (1973).
13. Effect of Ficoll-Hypaque separation on activation and DNA synthesis of human blood lymphocytes. Farnes, P., Barker, B.E. In: *Regulatory Mechanisms in Lymphocyte Activation*, (Ed. Lucas, D.O.) *Academic Press*, New York, 368-370 (1977).
14. Comparison of lymphocyte function after isolation by Ficoll-Hypaque flotation or elutriation. Berger, C.L., Edelson, R.L. *J Invest Dermatol*, **73**, 231-235 (1979).
15. Lack of subclass selection during density purification of human lymphocytes. Häyry, P., Tötterman, T.H., Ranki, A. *Clin Exp Immunol*, **28**, 341-346 (1977).
16. Special requirements for isolation of purified lymphocyte populations in infected patients. Dougherty, P.A., Balch, C.M. *Fed Proc*, **40**, 1120 (1981).
17. Changes in Ficoll-Hypaque gradients with advancing stage of lung cancer. Check, I.J., Hunter, R.L. *Fed Proc*, **38**, 1224 (1979).

18. Contamination of mononuclear cell suspensions obtained from cancer patients by the Bøyum method. Currie, G.A., Hedley, D.W., Nyholm, R.E., *et al. Brit J Cancer*, **38**, 555-556 (1978).
19. Non-lymphoid cells obtained by the Bøyum technique and their significance in cancer patients. Kluin-Nelemans, J.C., van Helden, H.P.T. *J Clin Lab Immunol*, **4**, 99-102 (1980).
20. Comparison of rates of virus isolation from leukocyte populations separated from blood by conventional and Ficoll-Paque/Macrodex methods. Howell, C.L., Miller, M.J., Martin, W.J. *J Clin Microbiol*, **10**, 533-537 (1979).
21. Gradient separation of normal and malignant cells. II. Application to *in vivo* tumour diagnosis. Minami, R., Yokota, S., Teplitz, R.L. *Acta Cytol*, **22**, 584-588 (1978).
22. A quick method for concentration and processing cancer cells from serous fluids and fine-needle nodule aspirates. Elequin, F.T., Muggia, F.M., Ghossein, N.A., *et al. Acta Cytol*, **21**, 596-599 (1977).
23. The comparative diagnostic accuracy of cancer-cell detection obtained with Ficoll-Hypaque gradient separation and standard centrifugation technics on body-cavity fluids. Katz, R.L., Lukeman, J.M. *Amer J Clin Pathol*, **74**, 18-24 (1980).
24. Enhancement of human amniotic cell growth by Ficoll-Paque gradient fractionation. Chang, H.-C., Jones, O.W., Bradshaw, C., *et al. In Vitro*. **17**, 81-90 (1981).
25. A simple method for the isolation of murine peripheral blood lymphocytes. Chi, D.S., Harris, N.S. *J Immunol Methods*, **19**, 169-172 (1978).
26. Isolation of various canine leucocytes and their characterization by surface marker analysis. Ho, C.K., Babiuk, L.A. *Immunol*, **35**, 733-740 (1978).
27. Lymphocyte isolation, rosette formation, and mitogen stimulation in rhesus monkeys. Taylor, D.W., Marchette, N.J., Siddiqui, W.A. *Develop Comp Immunol*, **2**, 539-546 (1978).
28. The bovine lymphoid system: Binding and stimulation of peripheral blood lymphocytes by lectins. Pearson, T.W., Roelants, G.E., Lundin, L.B., *et al. J Immunol Methods*, **26**, 271-282 (1979).
29. Acid -naphthyl acetate esterase: presence of activity in bovine and human T and B lymphocytes. Yang, T.J., Jantzen, P.A., Williams, L.F. *Immunol*, **38**, 85-93 (1979).
30. Humoral and formed elements of blood modulate the response of peripheral blood monocytes. I. Plasma and serum inhibit and platelets enhance monocyte adherence. Musson, R.A., Henson, P.M. *J Immunol*, **122**, 2026-2031 (1979).
31. Comparative study of six methods for lymphocyte isolation from several mammalian sources and determination of their carbohydrate composition. Hueso, P., Rocha, M. (Article in Spanish) *Rev Esp Fisiol*, **34**, 339-344 (1978).
32. Separation of porcine blood cells by means of Ficoll-Paque. Wittman, G. (Article in German) *Zbl Vet Med B*, **27**, 253-256 (1980).
33. A comparison of the methods used for the separation of fish lymphocytes. Blaxhall, P.C. *J Fish Biol*, **18**, 177-181 (1981).
34. Separation of human peripheral blood monocytes on continuous density gradients of Polyvinylpyrrolidone-coated silica gel (Percoll). Brandslund, I., Møller-Rasmussen, J., Fisker, D., *et al. J Immunol Methods*, **48**, 199-211 (1982).

35. Efficient separation of human T lymphocytes from venous blood using PVP-coated colloidal silica particles (Percoll). Feucht, H.E., Hadam, M.R., Frank, F., *et al.* *J Immunol Methods*, **38**, 43-51 (1980).
36. An improved technique for the isolation of lymphocytes from small volumes of peripheral mouse blood. Mizobe, F., Martial, E., Colby-Germinario, S., *et al.* *J Immunol Methods*, **48**, 269-279 (1982).
37. EC Guide to GMP (Good Manufacturing Practice), annex 1 "Manufacture of Sterile Medicinal Products.
38. United States Pharmacopeia. *Recommendations for ancillary materials*, chapter <1043>.

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Product		Quantity	Code no.
Ficoll-Paque PREMIUM	A density gradient medium manufactured in a GMP-compliant environment and certified to ISO 13485 standard; intended for the preparation of stem cells from blood and bone marrow for clinical applications.	6 × 100 ml	17-5442-02
		6 × 500 ml	17-5442-03
Ficoll-Paque PLUS	A density gradient medium for <i>in vitro</i> isolation of lymphocytes from whole human peripheral blood.	6 × 100 ml	17-1440-02
		6 × 500 ml	17-1440-03
Ficoll PM400	A hydrophilic polymer of high molecular weight for density gradient centrifugation.	100 g	17-0300-10
		500 g	17-0300-50
Percoll	A density gradient centrifugation medium for the isolation of cells, subcellular particles, and viruses.	250 ml	17-0891-02
		1 l	17-0891-01
Percoll PLUS	A silane-coated, silica-based medium for the isolation of various cell types for clinical research applications.	250 ml	17-5445-02
		1 l	17-5445-01
Density Marker Beads	Small colored beads of accurately known densities, for calibration of density gradients of Percoll and Percoll PLUS.	10 vials	17-0459-01
Protein A Sepharose 6MB	For separation of cells by affinity chromatography.	10 ml	17-0469-01
CNBr Activated Sepharose 6MB	For separation of cells by affinity.	15 g	17-0820-01
Cytodex™ 1	Microcarriers for cell culture.	25 g	17-0448-01
Cytodex 3	Microcarriers for cell culture.	10 g	17-0485-01

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