


Purification of Antibody Fragments 抗体片段的纯化



imagination at work

Outline 概况

Introduction 介绍

Applications 应用

- Fab
- Dab
- scFv

Summary 总结



GE imagination at work




imagination at work

Development of MAb market

单抗市场的发展

Estimated sales of biologics, 2007-2014



What comes beyond MAb's?


哪种产品将超越单抗？



Source: GE Healthcare estimations based on market analysis and company annual reports

Antibody fragment market

抗体片段市场




- Fab market already existing 抗体片段市场已经存在
- Shift towards smaller antibody fragments in early phases 早期转向更小的抗体片段市场

Source: Development trends for therapeutic antibody fragments, Nature Biotechnology Vol 27, No 4, April 2009


Antibodies and antibody fragments

抗体和抗体片段



Affinity ligands for antibody fragments


针对抗体片段的亲和配基



Two heavy chains
两条重链

Two light chains
两条轻链
Kappa or Lambda

Chromatography media for fragments 针对抗体片段的层析填料



Two heavy chains
两条重链

Two light chains 两条轻链
Kappa or Lambda

Potential and challenges with antibody fragments 抗体片段的机遇和挑战

Potential 机遇

Challenges 挑战

- **Improved kinetics** in solid tumors
实体瘤中提高药物动力学
- **Enables binding** of targets not accessible to Ab
可结合靶点而不影响抗体
- **Higher diversity** in specific binding structures
特殊结合位点上具有更高的可变性
- **Advantages in manufacturing** – yeast and bacterial systems can be used
工业生产上的优势——可使用大肠杆菌和酵母系统


- **Lower stability and solubility**
可溶性和稳定性差
- **Sometimes low target retention**
某些片段靶向性差
- **Varying expression levels**
表达水平不同
- **Lack of generic purification protocols**
缺少通用的纯化方案



Properties of Protein L 蛋白L的属性

Protein L binds: 结合

- to the variable region of most subtypes of Ab kappa light chain
结合抗体kappa轻链可变区
- a wider range of Ab classes than Protein A or G
比蛋白A或蛋白G有更广的抗体种类
- Ab fragments 抗体片段 ; Fabs, Dabs and scFv
- without interfering with antigen binding 对于抗原的结合没有干扰



Outline 概况

Introduction 介绍

Applications 应用

- Fab
- Dab
- scFv

Summary 总结




Purification of human Fab 纯化人源Fab

Sample: Fab kappa derived from IgG₁ in *E. coli* cell culture supernatant

样品：源于大肠杆菌细胞培养表达的IgG₁的kappa片段

- Concentration in feed: 1 mg/ml 上样浓度： 1 mg/ml
- pI: 8.5 等电点： 8.5
- Mw: 48 kDa 分子量： 48kDa




Capto™ L: Design of Experiments (DoE)

Capto™ L实验设计 (DoE)

DoE:

- Wash pH: 4.5-7.5
- Wash [NaCl]: 40-460 mM
- Elution pH: 2.9-3.1



Constants: 非变量


Load 15 mg/ml (70% of DBC), residence time 4 minutes, wash volume 7 CV
15mg/ml上样（动态载量的70%），保留时间4min，冲洗7倍柱体积

Responses: 响应值

Yield and *E. coli* protein (ECP) content
收率及大肠杆菌蛋白含量

Capto™ L DoE results: ECP


Capto™ L实验设计结果：大肠杆菌蛋白



Start: 300 000 ppm


Capto™ L DoE results: Yield

Capto™ L实验设计结果：收率



14
19/09/2012


Capto™ L: Sweet spot plot



Purification of mouse Fab

鼠源Fab的纯化

Column 层析柱: HiTrap™ Protein L, 1 ml
Sample 样品: 2 mg mouse polyclonal Fab containing kappa light chain
(Jackson ImmunoResearch laboratories) 2mg 含有kappa 轻链的鼠源多克隆Fab
Elution 洗脱: pH 7.4-2.3 in 10 column volumes
System 系统: ÄKTA™ avant 25



Binding buffer: PBS (10 mM Na-phosphate, 140 mM NaCl, 2.7 KCl, pH 7.4)
Wash buffer: 25 mM sodium citrate, 25 mM sodium phosphate, pH 7.4
Elution buffer: 25 mM sodium citrate, 25 mM sodium phosphate, pH 2.3

Flow rates: 1 ml/min (equilibration, wash and elution), 0.25 ml/min (sample load)



Outline概况

Introduction介绍

Applications应用

- Fab
- Dab
- scFv


Summary 总结



GE imagination at work

Purification of a human Dab 人源Dab的纯化 Determination of optimal elution pH 优化洗脱pH

Column 层析柱: HiTrap™ Protein L, 1 ml
Sample 样品: 1 ml extract of *E. coli* expressing a domain antibody (Dab) (D115)
1ml 大肠杆菌表达的Dab(D115)提取物
Elution 洗脱: pH 6.0-2.3 in 20 column volumes
System 系统: ÄKTA™ pure



Binding buffer: PBS (10 mM Na-phosphate, 140 mM NaCl, 2.7 KCl, pH 7.4)
Wash buffer: 50 mM citrate, pH 6.0
Elution buffer: 50 mM citrate, pH 2.3

Flow rates: 1 ml/min (equilibration, wash and elution), 0.5 ml/min (sample load)





imagination at work



imagination at work

Purification of a human Dab 人源Dab的纯化 Scaling up using step elution 使用步级洗脱来放大

Column 层析柱: HiTrap™ Protein L, 5 ml
 Sample 样品: 127.5 ml extract of *E. coli* expressing a domain antibody (Dab) (D115)
 Elution 洗脱: pH step gradient; pH 6.0 and pH 2.6
 System 系统: ÄKTA™ pure



Binding/wash 1 buffer: PBS (10 mM Na-phosphate, 140 mM NaCl, 2.7 KCl, pH 7.4)
 Wash 2 buffer: 50 mM citrate, pH 6.0
 Elution buffer: 50 mM citrate, pH 2.6

Flow rates: 5 ml/min (equilibration, wash and elution), 2.5 ml/min (sample load)



Outline 概况

Introduction 介绍

Applications 应用

- Fab
- Dab
- scFv

Summary 总结



Challenging customer application 客户应用的挑战

BIBITEC GmbH is a German CMO, specialized in the production of recombinant proteins and MAb for use in clinical trials up to phase III using mammalian cells.

BIBITEC GmbH是一家德国的CMO公司，专业生产用哺乳动物细胞培养、用于三期临床试验的重组蛋白和单克隆抗体。

In this application, the challenge was to **purify scFv fusion protein (57 kDa) from transgenic animal plasma**.

从转基因动物中纯化scFv融合蛋白（57kDa）的应用挑战

The sample had high content of albumin, IgG as well as many other plasma proteins

样品中含有大量的白蛋白、IgG及其它血浆蛋白

Dynamic binding capacity ($Q_{b, 10\%}$): 23 ma/ml

动态结合载量: 23 mg/ml





imagination at work



scFv, ~28 kDa Bis-scFv, bispecific, ~55 kDa

First screening: Promising results 首先筛选：理想结果

Column 层析柱: HiScreen™ Capto L, 4.7 ml
Sample 样品: Precipitated plasma containing scFv fusion protein
Elution 洗脱: pH step gradient; pH 6.0 and pH 2.0
System 系统: AKTA™



Lanes

1. Feed (1/10) 上样
2. Flowthrough (1/10) 流穿
3. Wash pH 8 (pH8 冲洗)
4. Wash pH 4 (pH4 冲洗)
5. Elution (1/4) 洗脱
6. Mw-marker 分子量Mark

Binding/wash 1 buffer: PBS (10 mM Na-phosphate, 140 mM NaCl, 2.7 KCl), 100 mM glycine pH 8.0
Wash 2 buffer: 50 mM glycine, 50 mM citrate, pH 4.0
Elution buffer: 50 mM glycine, 50 mM citrate, pH 2.0
Flow rate: 1 ml/min

Target protein 目标蛋白





imagination at work



imagination at work

pH gradient: Potential for optimization pH梯度：优化潜力

Medium: Capto™ L, 2.4 ml (column Ø=10 mm)
Sample: Precipitated plasma containing scFv fusion protein
Elution: pH 4.0-2.0 in 10 column volumes
System: ÄKTA™



Lanes
1. Flowthrough 流穿
2. Wash pH 8 (pH8冲洗)
3. Wash pH 4 (pH4冲洗)
4. Gradient peak 1 梯度峰1
5. Gradient peak 2 梯度峰2

Binding/wash 1 buffer: PBS (10 mM Na-phosphate, 140 mM NaCl, 2.7 KCl), 100 mM glycine pH 8.0
Wash 2 buffer: 50 mM glycine, 50 mM citrate, pH 4.0
Elution buffer: 50 mM glycine, 50 mM citrate, pH 2.0
Flow rate: 1 ml/min





imagination at work

Target protein 目标蛋白

Scale-up study 放大研究

Column 层析柱: Capto™ L, 34 ml (column Ø=26 mm)
Sample 样品: Precipitated plasma containing scFv fusion protein
Elution 洗脱: pH step gradient; pH 3.5 and pH 2.0
System 系统: ÄKTA™



Lanes
1. Feed 上样
2. Flowthrough 流穿
3. Wash pH 8 (pH8冲洗)
4. Wash pH 3.5 (pH3.5冲洗)
5. Elution 洗脱
6. Mw-marker 分子量Marker

Binding/wash 1 buffer: PBS (10 mM Na-phosphate, 140 mM NaCl, 2.7 KCl), 100 mM glycine pH 8.0
Wash 2 buffer: 50 mM glycine, 50 mM citrate, pH 3.5
Elution buffer: 50 mM glycine, 50 mM citrate, pH 2.0
Flow rate: 7.4 ml/min



imagination at work

Target protein 目标蛋白

Scale-up study放大研究

Contaminant levels	Plasma	Run 1	Run 2	Reduction
Albumin (ppm)	3.7×10^7	2.8×10^4	2.8×10^4	3 log
IgG (ppm)	1.8×10^8	7.6×10^3	7.9×10^3	4.5 log
Plasma proteins (ppm)	1.3×10^9	1.6×10^4	1.8×10^4	5 log
Pr L leakage (ppm)		≤ 2	≤ 2	

The yield was at least 89% for both runs
两次实验的收率至少达到89%



Conclusions: scFv application example 结论：scFv应用实例

- **Dynamic binding capacity: 23 mg/ml**
动态结合载量：23mg/ml
- **Challenging sample**
具有挑战性的样品
- **Effective wash step developed**
开发有效的清洗步骤
- **High yields obtained: >89%**
高回收率：>89%



Bis-scFv, bispecific, ~55 kDa



scFv, ~28 kDa



Outline概况

Introduction介绍

Applications应用

- Fab
- Dab
- scFv

Summary 总结



GE imagination at work

Summary总结

- Capto™ L is suitable for purification of antibody fragments containing kappa light chain;

Capto™ L适于含有kappa轻链的抗体片段的纯化

- Fab
- Dab
- ScFv

- Capto L purification show high yields and purity

Capto™ L纯化表现出较高的收率和纯度



imagination at work



imagination at work

Thank you!



ÄKTA, HiTrap, HiScreen, MabSelect, Capto and Sepharose are trademarks of GE Healthcare companies. GE and GE Monogram are trademarks of General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

© 2012 General Electric Company – All rights reserved.

GE Healthcare Bio-Sciences AB, a General Electric Company.

www.gelifesciences.com

GE Healthcare Bio-Sciences AB
Björkgatan 30
SE-751 84 Uppsala
Sweden



GE imagination at work



imagination at work