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
Cidra Puerto Rico

膜铸造与切割
Milliflex, Steritest 介质,
Dip tester, Petri dishes,
安瓿

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
Danvers

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Millipore Product Range: from a \$2.00 pipette tips to \$5 million custom-engineered systems

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



Why Integrity Test?

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Common Sense 通常理解

- Filtration is often the most critical step in an operation
- 过滤通常是操作的关键步骤
- Confirmation of manufacturers specifications
- 确认制造规格
- Detecting leaks due to o-rings, gaskets, seals
- 检测O形环, 垫圈, 密封垫的泄漏
- Assuring the correct pore size filter
- 确认正确的过滤孔径
- Assuring integrity before sterilization
- 确认灭菌前完整性
- Assuring integrity after steaming or autoclaving
- 确认蒸汽和消毒锅灭菌后完整性

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Guidance for Industry
Sterile Drug Products
Produced by Aseptic Processing —
Current Good Manufacturing Practice

Additional copies are available from:
Office of Training and Communication
Division of Drug Information, HFD-240
Center for Drug Evaluation and Research
Food and Drug Administration
5600 Fishers Lane
Rockville, MD 20857
(Tel) 301-827-4573
<http://www.fda.gov/cder/guidance/index.htm>

or
Office of Communication, Training and
Manufacturers Assistance, HFM-40
Center for Biologics Evaluation and Research
Food and Drug Administration
1401 Rockville Pike, Rockville, MD 20852-1448
<http://www.fda.gov/cber/guidelines.htm>
(Tel) Voice Information System at 800-835-4709 or 301-827-1800

U.S. Department of Health and Human Services
Food and Drug Administration
Center for Drug Evaluation and Research (CDER)
Center for Biologics Evaluation and Research (CBER)
Office of Regulatory Affairs (ORA)

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
→ Sterilization induced failures

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→ Stress induced failures

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→What Are The Integrity Testing Choices?

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
MILLIPORE

Regulatory & quality organizations need data from both to assure reliable and predictable filter performance
法规和质量管理需要数据确保可靠和预知过滤器性能

Destructive vs. Non-destructive Testing破坏性和非破坏性测试



There must be a correlation between a destructive and a non-destructive integrity test
破坏性和非破坏性测试两者必须有关联

- Regulatory requirement
- 法规需要
 - FDA Aseptic Guidelines
 - FDA无菌指南
- Validation justification for the use of a test
- 测试应用的验证理由

Integrity Testing Correlations

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Bubble point can have a direct correlation

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MILLIPORE

Diffusion & other tests can have an "go - no go" correlation

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
What Non-Destructive Integrity Tests are Available?

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
- Diffusion 扩散
- Forward Flow 前向流
- Pressure hold / decay 压力保持 / 衰减

Bubble Point Introduction 起泡点介绍

Bubble point is the pressure at which gas displaces liquid from the largest set of filter pores and flows rapidly through the filter

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
→ Bubble point indicates the magnitude of the forces holding liquid in the filter structure

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
→ The oldest non-destructive integrity test

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Bubble Point Test Applicability

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Useful for qualified manual testing of filters less than approximately 2000 cm² (~2 sq.ft.)

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Useful for automatic integrity testers for filters with gas flow rates of up to 100 ml/min or limit of the automatic testers qualifications

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Bubble Point Equation 起泡点等式

The bubble point is expressed as: 泡点可表达为

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



Must know the effects of the changes, perform flushing or conduct product bubble point qualification

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


Diffusion Testing Introduction

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- Gas dissolves in liquid held in the pores of a fully wetted membrane filter. 气体溶解在完全湿润的滤器膜孔内的液体中
- A pressure differential will give a different gas concentration across the filter. 膜压差造成跨膜气体浓度差
- Results in gas flow through the liquid dissolved in the filter pores. 结果使膜孔溶解的气体流出

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Diffusion Testing Equation 扩散测试方程

$$\text{Diffusion} = \frac{K \cdot (P_1 - P_2) \cdot A \cdot \rho}{L}$$

where:


K = Diffusivity / Solubility
coefficient 扩散因子

P₁, P₂ = Pressure difference across
the system 系统两边的压力

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
Bubble Point detection becomes harder as membrane area increases

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Need as much information about the filter as possible

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- Use diffusion testing in combination with bubble point testing 使用起泡点与扩散测试相结合

Diffusion Test Applicability

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For filters with gas flowrate of >100 ml/min

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- filter type 过滤器类型
- wetting liquid 湿润液体
- test gas 测试气体
- Temperature 温度

■ maximum acceptable flow rate 最大接受流速

Factors Affecting Diffusion Tests

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Wetting湿润度

- Lower or inconsistent film thickness increases diffusion
- 低的和不均匀的薄厚会增加扩散

Fluid流体

- Organic solvents have high diffusion flowrates
- 有机溶剂有较高的扩散流速

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
- Carbon dioxide has higher gas flow than air or nitrogen
 - 二氧化碳比空气或氮气有较高的扩散流速

Temperature温度

- 80% higher flowrate at 60 degrees
- 在 60 摄氏度流速高80%

MILLIPORE

“Diffusion Testing” Summary

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Diffusion testing measures gas flow 扩散测试测量气体流速

At bubble point, gas flows rapidly through the largest pore(s)

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
Downstream gas flow (or upstream pressure loss) can come from both diffusion through the wetted pores PLUS the gas flow through the “dry” pores which have reached their bubble point

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- **Alcohol / Water mixtures are often used due to:**
- 酒精/水的混合物使用经常会导致
 - Cost成本高
 - Flammability有燃烧的危险
 - high diffusion高的扩散流
- **Secure areas and equipment needed**需要可靠的区域及设备
- **Can be difficult to perform in-situ**很难在线操作
 - **Some concerns with residual test solution remaining in filter holder / pipework after testing**
 - 在测试后,测试液体会残存在过滤器/管路中
 - **Residual solution should be removed before filter sterilization**
 - 灭菌前需要去除残存的液体


HydroCorr Testing

Hydrocorr测试


- Force water under pressure onto the surface of the filter.
- 把水灌在膜表面
- Applied pressure is $<$ integral filter intrusion pressure
- 加压至小于膜的水浸入压
- The water pressure is maintained but not increased.
- 压力维持但不再加压
- Pleated structure compacts
- 折叠式结构进一步压缩
- If integral, the upstream pressure will not drop below a preset value.
膜如果是完整的,上游的压力不会太快下降

MILLIPORE


HydroCorr Test

Resistance to water intrusion 水浸入法测试

HydroCorr Test Results测试结果

- After initial stabilization, “apparent” water flow for an integral filter is very low
- 在最初的稳定后,对于完整膜来说,水流是非常小的
- Specifications are correlated to bacterial challenge
- 水流指标是与细菌挑战有关
- For integral filter there is no “actual” water flow
- 实际上完整膜,没有实际的水流
- Filter remains dry and does not require drying between testing and use
- 在测试与使用中,过滤膜维持干燥也不需要干燥

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

HydroCorr Test Considerations 测试的注意点

- Must have a leak-free system 必须在一密闭系统内进行
- Filters must be dry and clean 过滤芯必须干燥且干净
 - It is relatively easy to Hydrocorr test “very hydrophobic” membranes such as PTFE
 - 对于疏水性很强的滤芯,例PTFE材质的膜,水浸入法是相当容易
 - Successfully testing “less hydrophobic” membranes such as PVDF requires ideal conditions that can be difficult to achieve in actual biopharmaceutical processes
 - 在实际的生物制药工艺中,测试一些疏水性不强的滤膜,例PVDF膜是比较困难的
- Temperature changes should be minimized
- 尽量减少温度的变化
- Reference testing is important
- 相关的参考比较重要
- Must have a complete validation package
- 必须有一完整的验证文本
- Values must be comprehensively correlated to microorganism challenge
- 数据必须与细菌挑战相关

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

Summary总结

- Bubble point and Diffusion (or “Forward flow”) are BOTH valid integrity tests as per regulatory agencies.

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- Bubble point provides a direct correlation to bacterial retention, the critical performance characteristic

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Automated Integrity Testers

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MILLIPORE

Why Do People Use Automatic Testers

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
- **Revision of annex 1 to EC Guide to GMP for sterile medicinal products (1997)**
- 欧盟对于无菌药品GMP指南的附件一(1997版)
 - “The integrity of the sterilized filter should be verified before use and should be confirmed immediately after use by an appropriate method such as a bubble point, diffusive flow or pressure hold test”
 - 除菌过滤膜应该在使用前及使用后立即采取合适的方法确认其完整性,可以采用泡点,扩散流或压力保持的方法.
- **FDA Guideline on Sterile Drug Products Produced by Aseptic Processing (1987)**
- FDA对于无菌操作生产的无菌药品的要求(1987版)
 - “Normally, integrity testing of the filter is performed after the filter unit is assembled and sterilized prior to use.”
 - 完整性测试通常在过滤器安装,灭菌后使用前进行
- **PDA Technical Report # 26**
- **PDA技术报告26版**
 - “It generally is regarded as a cGMP requirement that filters or filter systems routinely be integrity tested both prior to and after use.”
 - 现行的GMP要求过滤器及过滤系统在使用前及使用后均需要进行完整性测试

Automated Integrity Testing

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MILLIPORE

A. During integrity testing gas molecules leave the system due to diffusion or bulk flow
在完整性测试中,气体流失的原因是扩散及穿孔后的流失

B. Fewer gas molecules upstream results in lower pressure.

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C. All else being constant, the change in pressure can be used to determine the gas flow rate

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- Key Point - The pressure in a closed system is directly related to the number of gas molecules in the system.
- 关键点:在密闭系统内,压力直接与气体的分子数有关

Relating Pressure Decay to Gas Flow Rate

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P_{atm} = atmospheric pressure 大气压

V_u = housing volume (cc) 外壳体积

How Automatic Testers Determine Integrity

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
- Dependent on upstream volume
- 依据上游体积
- Dependent on gas flow curve properties of the filter
- 依据过滤器的气体流量曲线
 - Diffusion is gas flow measurement at a single pressure
 - 扩散是在单一压力下测量的气体流
 - Bubble point is determined by measuring gas flow multiple times at increasing pressures
 - 泡点是在不断加压的情况下多次测量气体流的情况下而得的。

Upstream Volume Measurement

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
→ $P1V1 = P2V2$

- Internal tank is V1
- 内部储罐是V1
- P1 is measured
- P1是测量而得
- Valve between internal tank and housing opens
- 在内部储罐及外壳之间的阀门是开的
- P2 is measured
- P2是测量而得
- V2 is calculated
- V2是计算而得

Gas Flow Curve 气体流量曲线


- Diffusion and Bubble point together make up a "gas flow curve"
- 扩散及泡点共同组成了一个“气体流量曲线”
 - If flow is due to diffusion, it will increase linearly vs. delta P
 - 如果气体流量是由扩散而造成,流速与压降成正比
 - If flow is a combination of diffusion and bulk flow, it will increase exponentially vs. delta P
 - 如果气体流量由扩散流及穿孔流共同组成,流速与压差成指数关系

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Diffusion / Bubble Point Procedure

Hydrophilic 0.2 micron filter

Effect of Temperature on Automated Testing

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- Temperature CHANGES during a pressure decay test will have MAJOR impact on test results.
- 在压力衰减测试中,温度变化对结果的影响是非常大的.
 - Ideal Gas Law - $PV = nRT$ 理想气体定律
 - As Temperature increases Pressure increases
 - 温度上升时压力上升
 - As Temperature decreases Pressure decreases
 - 温度下降时压力下降
- If the temperature changes even 1°C . during the test, the pressure change WILL NOT be solely due to gas flow
- 在测试中,温度变化哪怕仅1度,压力的变化将不在仅归于气流
 - Additional pressure loss or pressure gain are test error
 - 额外的压力损失或增加会发生测试错误

Temperature Changes温度变化因素

- Temperature changes can be due to:
- 温度变化是由以下因素引起
 - Inadequate stabilization after steaming
 - 灭菌后,冷却时间不够
 - Touching the housing碰到了外壳
 - Heating the automated tester加热了自动完整性测试仪
 - Environmental changes during testing在测试中,环境温度变化
 - Expansion and contraction膨胀及压缩

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- Important to look at the whole printout, not just the words “PASSED” or FAILED”
- 浏览整个打印报告是非常重要的,而非仅看结果”通过”或”失败”
- Check for:核实
 - Changes in gas flow rate 气体流速的改变
 - Temperature changes during the test 测试中温度的变化
 - Very low gas flow rates 气体流速非常低
 - Increasing temperature 温度的提高
 - Improperly closed valves 不正确的关闭阀门
 - Marginally high gas flow rates 非常高的气体流速
 - Leaks 泄露
 - Product specifications 产品的规格

Summary总结


- Main advantage of automatic integrity testers is ability to easily test sterilized systems
- 自动完整性测试仪的优点是能够容易的测试灭菌后的系统
- Upstream pressure decay measurements are related to actual downstream gas flow
- 上游压力衰减与下游气体的流速相关
- Automated integrity test systems are VERY sensitive to temperature changes
- 自动完整性测试仪对于温度的变化非常敏感
- Interpreting the test result is more than just looking for “PASSED” or “FAILED”
- 看结果时不应只看”通过”或”失败”

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



Administrator

Administrator w/ Service privilege

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- US GMP 21 CFR Parts 210 & 211 美国GMP
- Appropriate written procedures...shall be established and followed. Such procedures shall include **validation of any sterilization process**
- 适当的书面规程...应当被建立和遵循。这些规程应当包括任何无菌工艺的验证。
- EU GMP Annex 1 Sterile Medicinal Products 欧盟 GMP
- All **sterilization processes** should be **validated**
- 所有的无菌工艺都应当被验证
- Australian TGA GMP 澳大利亚TGA GMP
- **Filtration processes** used as the sterilizing step for products should be **validated**
- 作为产品除菌步骤的过滤工艺应当被验证
- Health Canada GMP 加拿大GMP
- Documented evidence is available establishing **validation** and validity **of each sterilization process**
- 每个无菌工艺的验证及其有效性都应当建立书面化的证据

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- **FDA Aseptic Processing Guidelines (1987 replaced by new 2004 - www.fda.gov)**
- **FDA无菌工艺指南（2004新版）**
- **Correlate filter performance with filter integrity testing**
- **Include microbiological challenges to simulate ‘worst case’ production conditions**
- 过滤器性能与完整性测试的关联应当包括微生物挑战以模拟“最差”的工艺条件”。
- **PDA Technical Report 26 (1998, www.pda.org) PDA技术报告26**
- **“Early, careful screening of potential filter types and configuration can result in fewer technical and regulatory problems, fewer delays, more efficient product processing, and greater sterility assurance”**
- 尽早的，仔细的选择过滤器类型及配置可以有效的减少技术及法规方面的问题，减少延迟，使生产工艺更为有效，提供更高的无菌保障。
- **ISO/DIS 13408-2 Aseptic Processing (2003, www.iso.ch) ISO/DIS 13408-2无菌工艺**
- **Filter Pre-selection shall take into account chemical and physical characteristics of the filter, as established by the filter manufacturer.**
- **Bacterial retention performance of filters shall be validated in a fluid-specific manner or for fluid groups under worst case conditions in production.**
- 过滤器的预选应当考虑到过滤器生产商所确定的滤器的化学物理特性。过滤器的细菌截流特性应当通过生产中特定的料液或者最差条件下的料液进行验证。

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




MILLIPORE

Contents

Product Summary	7
1.1 Introduction	7
1.2 Product Description	7
1.3 Catalogue Number	8
1.4 Product Specifications	8
1.5 Certification	9
1.6 Test Summary	11
1.7 Chemical Compatibility	13
Millipore Qualification Test Results	19
2.1 Organism Retention	19
2.1.1 Bacterial Retention/Integrity After 1 Steam-In-Place Cycle Of Thirty Minutes, 135°C (31 psig or 2.14 bar gauge steam)	19
2.1.2 Bacterial Retention/Integrity After 30 Steam-In-Place Cycles Of Thirty Minutes, 135°C (31 psig or 2.14 bar gauge steam)	21
2.2 Integrity Test Values	23
2.2.1 Diffusion, Bubble Point, and Reverse Integrity	23
2.3 USP Biological Safety	25
2.3.1 USP Class VI Biological Tests For Plastics	25
2.3.2 USP General (Mouse) Safety Test	32
2.3.3 USP Bacterial Endotoxin	37
2.4 Water Flow Rate	42
2.5 Hydraulic Stress Resistance	47
2.5.1 Hydraulic Stress at 25°C	47
2.5.2 Hydraulic Stress at 80°C	49
2.5.3 Hydraulic Stress at 135°C	50
2.6 Thermal Stress Resistance	52
2.6.1 30 Steam-In-Place Cycles at 135°C	52
2.6.2 30 Autoclave Cycles at 126°C	53
2.7 Cleanliness	53
2.7.1 Gravimetric Extractables in Water	54
2.7.2 USP Oxidizable Substances	55
2.7.3 Particle Releasing	56
2.7.4 Fiber Shedding	58
Millipore Qualification Test Methods	63
3.1 Organism Retention	63
3.1.1 Bacterial Retention/Integrity After 1 Steam-In-Place Cycle Of Thirty Minutes, 135°C	63
3.1.2 Bacterial Retention/Integrity After 30 Steam-In-Place Cycles Of Thirty Minutes, 135°C	63
3.2 Integrity Test Values	65
3.2.1 Diffusion, Bubble Point, and Reverse Integrity	65
3.2.2 Cartridge Wetting Procedure	67
3.3 USP Biological Safety	68
3.4 Water Flow Rate	69
3.5 Hydraulic Stress Testing	71
3.5.1 Hydraulic Stress Resistance at 25°C	71
3.5.2 Hydraulic Stress Resistance at 80°C	72
3.5.3 Hydraulic Stress Resistance at 135°C	72
3.6 Thermal Stress Resistance	73
3.6.1 30 Steam-In-Place Cycles at 135°C for 30 minutes	73
3.6.2 30 Autoclave Cycles at 126°C for 60 minutes	74
3.7 Cleanliness	75
3.7.1 Gravimetric Extractables in Water	75
3.7.2 USP Oxidizable Substances	76
3.7.3 Particle Releasing	77
3.7.4 Fiber Shedding	78

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


FIG. 1 Test Set-Up for Bacteria Retention Testing

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
Cartridge Filters —
0.1 μm Hydrophilic Durapore Membrane

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