ÄKTA™ avant

Inspired to Accelerate your process development
有灵感的加速您的工艺开发

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 imagination at work

Designed for process development and scale-up. 工艺开发与中试放大

Target discovery, HTS → Drug substance development → Toxicity, pre-clinical studies → Clinical Phase I → Clinical Phase II → Clinical Phase III → Approval Phase → Post-marketing studies

First external advice → Confidence in Safety → Confidence in Mechanism → 新药证书 → 正式生产文号

ÄKTAavant Targeting

Process Development and Scale-up phase

ÄKTAavant is based on customers' need for **speed, security and scalability** during process development and Scale-Up.

Workflow

- Expression
- Clarification
- Purification
- Formulation
- Filling

 imagination at work

Two versions for maximized performance 两款性能最佳的系统

ÄKTA™ avant 25

- Screening in small columns
小柱子进行筛选实验
- Method optimization
方法优化
- Challenging purifications
纯化挑战性蛋白
- Small amounts of sample
小量的样品
- High res. prep. Purification
高分辨率制备纯化
- Optimized for column size i.d. 5-16 mm
在直径5-16mm柱上进行优化

ÄKTA avant 150

- Scale-up process
工艺放大
- Robustness testing
稳健性实验
- Trouble shooting (scale-down)
问题解决 (缩小规模)
- Fine-tuning of process
工艺的精细调节
- Optimized for column size i.d. 10-50 mm
在直径10-50mm柱上进行优化

A complete workflow for Process development 一个工艺开发的完整流程

Predict 筛选 Optimize 优化&Scale-up放大 Produce 生产

Predict 筛选
PreDicator™
HTPD
- Media screening
- 填料筛选
- Condition screening
- 条件筛选

Optimize 优化&Scale-up放大
ÄKTA™ avant 25
Screening 筛选
Method optimization
方法优化
Column size i.d. 5-26 mm

ÄKTA avant 150
Scale-up process 放大工艺
Fine-tuning 精细调节
Robustness testing
稳定性测试
Column size i.d. 10-50 mm

Produce 生产
ÄKTApilot™ ÄKTA ready
Scale-up optimized process
最大化的工艺
Sanitary design
卫生设计
Column sizes from i.d. 50 mm

 imagination at work

5 /
GE 28-9577-83 /
Version 4.0

ÄKTA™ avant Setting the future industrial standard ÄKTA™ avant, 设定新的行业标准

 Inspired speed 卓越速度

 imagination at work

GE 28-9577-83 /
Version 4.0


Speed 速度

ÅKTA™ avant has the capacity to significantly shorten the process development time, allowing you to get your products to the market faster
AKTA avant 可以显著性地缩短您工艺开发的时间，加快您的新产品上市

Maximum information, minimum experiments with DoE

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Automation with BufferPro
BufferPro缓冲液自动配制提高自动化程度

Ease-of-use with Method editor
易用的方法编辑

7 /
GE 28-9577-83 /
Version 4.0

What is Design of Experiments (DoE)? 什么是实验设计?

DoE is a systematic way of changing process inputs and analyzing the resulting process outputs in order to quantify the cause and effect relationship between them while using a minimum number of runs.

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
Basis for setting acceptable and realistic variability ranges for control space
将参数控制区域设定在可接受和可实现的变化范围

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



DoE arranged approach:
Vary all simultaneously

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
Example - No interaction
Effect of pH independent of T

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

Example - Strong interaction
Effect of pH reversed at high T

$$y = \beta_0 + \beta_1 x_1 + \beta_2 x_2 + \beta_{12} x_1 x_2$$

12 /
GE 28-9577-83 /
Version 4.0

The Model

X : known input - experimental settings, factors

Y : measured - response

b : unknown constant - determined by the regression model

Constant Term = response of Y when factors are 0

$$Y = b_0 + b_1X_1 + b_2X_2 + b_3X_3 + b_{12}X_1X_2 + b_{13}X_1X_3 + b_{23}X_2X_3 + b_{11}X_1^2 + b_{22}X_2^2 + b_{33}X_3^2$$

Linear Terms = main effects

Fractional factorial design
Screening or Robustness testing

Two-way Interaction Terms

Full or Fractional factorial design
Screening

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


DoE Screening 筛选

DoE Optimization 优化

Operation and DoE
Robustness 运行和稳健性

Or in QbD terms...(FDA guidelines)...

Systematic approach to development:

- Systematic approach to development: 系统性的工艺开发方法
- Begins with predefined objectives 起始设定好目标
- Emphasizes product and process 强调产品和工艺
- Understanding and process control 理解和过程工艺控制
- Based on sound science and quality risk management 以科学和质量风险管理为基础

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Product 产品

Acceptable variability in key and critical quality attributes 关键质量影响因素可接受的变化范围

Process characterization studies 工艺参数研究

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
Mw ~11 000

pI ~ 5.2

8M Urea is needed to solubilize inclusion bodies of r-Pro-insulin 需要用8M的尿素来溶解重组胰岛素原包涵体

Cation, Anion exchangers and Multimodal resins may be used 阳离子、阴离子交换填料和多模式配基填料供选择

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ÄKTA™ avant 25 & ÄKTA avant 150

-The possibility to automate process development

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


18 /
GE 28-9577-83 /
Version 4.0

Aim of study 研究的目标

Fast and efficient development of a r-Pro-Insulin capture step through High Throughput Process Development (HTPD) to scale up

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


Sample: 200 µl of total protein conc. 4.5 mg/ml

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i.e. 4 x diluted crude sample to adjust conditions

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
Result: No binding of r-Pro-Insulin to AIEX and Capto™ adhere

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
Conclusion: Capto MMC at 150 mM salt is chosen

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Conclusion:
Highest binding capacity = ~16mg/ml media at pH ~5.2 & 0-150 mM NaCl assuming that r-Pro-Insulin concentration is ~30% of total protein conc.

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


Conclusion: 1000 mM salt is not needed for elution

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
Experimental:

Column: Tricorn™ 5/50 - Capto™ MMC
Res time: 5 min (0.2 ml/ml).
Sample: Adjusted to pH 5.2 with HAC. 2.5 C.V. crude sample was loaded
Equilibration: 50 mM Na acetate + 150 mM NaCl pH 5.2*
Step elution: Study pH 6.2-8 and 150-750 mM NaCl using Acetate/Phosphate/Formate bufferPro (CIEX 0-1M NaCl pH2-7)*

No. of experiments: 11 (3 center point included)
实验数目：11 (包括3个中心点重复实验)


*Buffers in 8 M urea (缓冲液含8M的尿素)

30 /
GE 28-9577-83 /
Version 4.0

Result: r-Pro-Insulin yield as a function of pH & NaCl concentration

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imagination at work

31 /
GE 28-9577-83 /
Version 4.0

Suggested loading and elution conditions for capture step of r-Pro-Insulin

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Resin: Capto™ MMC
Loading: Adjust pH of crude sample to 5.2
Load 2.5 x Column Volumes
Step elution: 50 mM Phosphate buffer pH 8.0
150 mM NaCl
8 M Urea

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
No. of experiments 实验数目: 11 (3 center point included, 包括3个中心点)

Conclusions: Stable results supports that these conditions will give robust purity results

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
Design of Experiments allows the maximum amount of information to be obtained from a minimum number of experiments. By improving efficiency, both time and money are saved. 实验设计通过最小的实验数量获得最大的实验信息量。通过提高效率节省金钱和时间。

- Design of Experiments (DoE) is integrated in UNICORN™ 6

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- Reach the same result in 16-32 experiment: as you would with 128 experiments using the traditional approach!

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- Completely automated, advanced on-line buffer preparation

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- Buffers are made in analogy with manual preparation: corresponding acids or bases used for pH adjustment

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- Actual mixing ratios and buffer recipes are provided

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
- Allows for easy and quick pH-scouting

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UNICORN 6 includes a new Method Editor, replacing the method wizard. No new text programming is required, new components immediately ready - save a day's work!

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- Contains all instructions used for controlling the chromatographic run

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- New user interface for easy viewing and editing of run properties

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- Use pre-defined methods / phases or program your own

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- Drag and drop from the phase library

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imagination at work

36 /
GE 28-9577-83 /
Version 4.0

System & sample pumps set up to save time

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ÄKTA™ avant sample and system pumps have the same technical design and run independently, no re-plumbing is necessary, - save a day's work per week!

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GE 28-9577-83 / Version 4.0

ÄKTA™ avant

Setting the future industrial standard

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Inspired scalability
卓越的放大性

GE 28-9577-83 / Version 4.0

Scalability 放大性

ÄKTA™ avant ensures smooth transition from product discovery to production through direct, reliable scalability using the same media

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Exploit the advantages of modern BioProcess™ media

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UNICORN™ 6 - A higher standard of intelligence.

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GE 28-9577-83 / Version 4.0

Increase your productivity with high flow media

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ÄKTA™ avant has flow rate and pressure specifications that support modern, high flow BioProcess™ media such as MobSelect™ and Capto™, allowing for smooth scale-up

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- Direct scalability to process scale, using the same media 用同样的填料直接放大到工业规模
- Increased productivity with high dynamic binding capacity at high flow 在高流速下具有高的结合容量增加了生产力
- Increased yield with rapid mass transfer 具有快速的传质增加了收率
- Reduced process time with high volume throughput 更大体积的流速可以减少工艺操作时间
- Cost-effective processing 工艺的成本效率

Full support of modern media...
Screening, column id 5-16 mm: 20-25 ml/min
Scale-up, columns id 16-50 mm: 120-130 ml/min

GE 28-9577-83 / Version 4.0

Any application without re-plumbing

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
ÄKTA avant allows higher flow rates and lower back pressure, supporting modern BioProcess media

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
GE 28-9577-83 / Version 4.0

GE 28-9577-83 / Version 4.0

Dedicated column formats for method screening 采用专用的柱子进行方法的筛选

Use HiScreen™ pre-packed columns - save a day's work and make scale up easy!

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GE 28-9577-83 / Version 4.0

Reliable, reproducible and easy to use!

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Pressure stable HiScale™ columns are compatible with modern BioProcess media

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GE 28-9577-83 / Version 4.0

Robust, intuitive design and easy to use!

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GE 28-9577-83 / Version 4.0

ÅKTA™ avant

Setting the future industrial standard

Inspired security

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GE 28-9577-83 / Version 4.0

GE 28-9577-83 / Version 4.0

Security 安全性

ÅKTA™ avant incorporates cutting-edge functionality critical to achieving real security throughout the development process

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
Cooled fraction collector for true sample integrity

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Reliable column protection with differential pressure measurement

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UniTag scanning for complete traceability

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53 /
GE 28-9577-83 /
Version 4.0

Fraction collector for optimal security 半导体致冷的收集器提高样品安全性

Hassle-free overnight running, outstanding reproducibility and sample security!

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- Temperature control (6 to 20°C) allows samples to be cooled and protected

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- Cassettes and rack trays available for micro titer plates up to 250 ml bottles

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- Load up to six cassettes in any combination

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- Automatic detection of tube rack type

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
53 /
GE 28-9577-83 /
Version 4.0

Precolumn pressure and delta pressure (ΔP) 柱前压和柱前柱后压差(ΔP)

For increased operational safety, the pressure difference over the packed media bed (ΔP) is continuously measured. 为了增加操作安全性, AKTA Avant 持续监测加在柱床上的压力降

- Pressure sensors integrated into the column valve 柱阀上整合了2个压力传感器
 - Before the column to protect column hardware 柱前的压力设置保护柱子
 - After the column, the pressure difference is calculated over the packed media bed (ΔP) 柱前柱后压力降可以直接计算出来用于保护柱床凝胶
- Automatic pressure-flow control to prevent overpressure 自动的压力流速控制模式来避免因堵柱而超压运行
- Pressure sensors are also connected to the system and sample pumps 系统泵和样品泵也设计了压力传感器

GE imagination at work

52 /
GE 28-9577-83 /
Version 4.0

Pressure limit – Definitions 压力限制-定义

Column hardware pressure limit is the pressure where the column breaks (normally = starts leaking). 柱子硬件的最大耐压是指引起柱子开始破裂的压力 (通常等同柱子开始渗漏压力)

p1

Maximal pressure over the packed bed, Δp is the pressure when the bed starts to compress. 柱床最大耐压, Δp 是指柱床开始压缩时候的压力

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
Pressure generated after the column 柱后产生的压力


The pressure at this position affects the column hardware pressure limit (= pre-column pressure) 此处压力影响柱子硬件最大耐压 (等同柱前压)

Pressure over the packed bed = delta column pressure, Δp 柱床压力=柱前后的压差 Δp

Pressure sensor 压力感应器

GE imagination at work

54 /
GE 28-9577-83 /
Version 4.0

Column handling optimized for sample security 柱子管理功能最大化地保护样品的安全

Remove the risk of failed runs due to column performance, lifetime expiry &/or misuse, by keeping track of the columns with the Column Logbook and UniTag labels

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Easy identification, automatic traceability and outstanding security of operation!

GE 28-9577-83 / Version 4.0

UniTags

GE 28-9577-83 / Version 4.0

ÄKTA™ avant

Setting the future industrial standard

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Buffer tray for convenient, large capacity storage.

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One wet side with modular set-up limits service down-time

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Column valve with built-in by-pass and flow direction change

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pH valve with built-in positions for pH calibration and storage. pH valve integrates calibration and storage functions

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Cooled, built-in fraction collector prevents local heating of samples.

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Swivel foot for easy access to all parts of the instrument. 旋转底座, 可锁定和解锁

Automatic sample application without re-plumbing. 自动上样的装置

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High precision system pumps with integrated air sensors.

GE 28-9577-83 / Version 4.0

System flow-path

GE 28-9577-83 / Version 4.0

Innovative setup and components 创新的设计和组件

Minimized down time during service with easily accessible front mounted units
易于接触正面部分组件, 可节省维护时间

- One wet working area easily accessed by simply by rotating the unit
和液体接触的工作区域通过旋转很容易接触到
- Components have same geometry, the modules can be quickly exchanged
同样结构的部件, 模块化很容易更换

Cabinet with easy access 	Instrument Display <ul style="list-style-type: none"> 5.7" TFT touch screen display Shows selected run data Has pause/continue functionality
UV Monitor UV <ul style="list-style-type: none"> No optical fibers accessible eliminating risk of damage - increased safety Easy to change flow cells: 2 and 10 mm True optical path length identified automatically 	Mixer <ul style="list-style-type: none"> Mixer sizes 0.6, 1.4 and 5 ml Integrated in-line filter
Optional Components <ul style="list-style-type: none"> Extra air sensor before sample pump or before column Up to 3 extra valves 	

Multi-functional units

61 /
GE 28-9577-63 /
Version 4.0

Sophisticated valve design 精巧的阀门设计

Completely redesigned valves - versatile system setup for easy method development
全新的阀门设计-全能的系统设计使工艺的开发更简单

- Increased productivity and yield is obtained with the column valve: it connects up to 5 different columns and has built-in reversed flow functionality
多柱包柱阀增加了生产力和产量, 新设计柱阀具有5个柱位, 并具备正反流控制功能
- Valve w. integrated pH-electrode enables easy calibration and automatic by-pass of electrode
pH 计阀门可以自动控制在线或离线, 并简化了校正工作

Injection Valve <ul style="list-style-type: none"> Multiple sample application techniques without re-plumbing System pump and Sample pump can be run simultaneously 	Column Valves <ul style="list-style-type: none"> Up to five column positions Built-in bypass position Built-in flow direction change (reversed flow) Built-in pressure sensors before and after column
pH Valve (5 Monitor) <ul style="list-style-type: none"> Built-in bypass of pH flow cell Built-in positions for pH calibration and storage Built-in bypass of flow restrictor 	
Inlet Valves <ul style="list-style-type: none"> Built-in air sensor useable with all inlets 	
Outlet Valve <ul style="list-style-type: none"> Waste and Frac positions Up to 10 outlet positions 	

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while giving you outstanding scalability and sample security 充分提供样品的安全性和可放大性

Thank you

GE Healthcare Bio-Sciences AB, a General Electric company.

GE Healthcare Bio-Sciences AB
Björkgräns 30
751 84 Uppsala
Sweden

PhaDicer, HiScreen, HiScale, HiTrap, Tricorn, ÄktaChrom, Capto, Sepharose, Mono Q, Äkta and ÄktaExplorer are trademarks of GE Healthcare Companies.
GE Imagination at work and GE Monogram are trademarks of General Electric Company.

All third party trademarks are the property of their respective owners.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

© 2010 General Electric Company - All rights reserved.
First presented October 2010

GE Imagination at work