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The purpose of the Design Qualification(DQ) is to provide documented evidence to prove the design of the M2 Downstream Process System and Automation System Project of Innovent Biologics(Suzhou) Co., Ltd. (hereafter referred to as "Innovent") complies with the requirements of the user requirement specification and relevant regulations.

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This DQ defines the purpose, procedure and acceptance criteria of tests to be executed for the mechanical, electrical, instrument and automation parts in M2 Downstream Process System and Automation System Project, which is supplied by SMP and delivered to Innovent.

编号 No.	参考 URS Ref. URS	子系统/模块 Subsystem/submodule	设计范围 Design scope
1	M2 填料匀浆罐用户需求标准（文件编号） M2 URS About Homogenize Tank System	纯化模块 Purification Module	工艺、容器、工艺管道、结构、电气、仪表 Process, vessel, process piping, structure, electrical, instrumentation
2	M2 纯化模块用户需求标准 M2 URS About Purification Module M2 原液分装系统用户需求标准 M2 URS About DS Preparation System	原液分装模块 DS Preparation Module	

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5 参考 REFERENCE

适用标准和法规	Applicable Standards And Regulations
(NMPA) 中国药品生产质量管理规范 (2010 年修订版)及其附录	(NMPA) Good Manufacturing Practice, revised in 2010, and related annex
欧盟药品法规第 4 卷 cGMP (人用和兽用药品) 及其附录	European Commission. The Rules Governing Medicinal Product in the European Union, Vol. IV, Good Manufacturing Practices. Medicinal products for human and veterinary use, and related annex
(FDA) 联邦法规第 21 篇第 210/211 部分, 成品药的现行生产质量管理规范	(FDA) 21 code of Federal Regulations Parts 210/211, Current Good Manufacturing Practice for Finished Pharmaceuticals
(FDA) 联邦法规第 21 篇第 11 部分, 电子记录和电子签名	(FDA) 21 code of Federal Regulations Parts 11, Electrical Record and Electrical Signature
(FDA)联邦法规第 21 篇第 177 部分子部 C177.2600-用于反复使用的橡胶制品	(FDA) 21 code of Federal Regulations Parts 11 subpart C section 177.2600-Rubber articles intended for repeated use
固定式压力容器安全技术监察规程 TSG R21-2016	Supervision Regulation on Safety Technology for Stationary P.V. TSG R21-2016
压力容器 GB/T150.1~150.4-2011	Pressure Vessel GB/T 150.1~150.4-2011
压力容器焊接规程 NB/T 47015-2011	Welding Specification for Pressure Vessel NB/T 47015-2011
承压设备无损检测 NB/T 47013-2015	Nondestructive Testing of Pressure Equipment NB/T 47013-2015
机械安全机械电气设备第一部分: 通用技术条件 GB5226.1-2008	Safety of machine---Electrical equipment of machines---Part 1:General requirements GB 5226.1-2008
机械安全 防护要求 GB12265.2-2000	Safety of machine Safety distance to prevent danger zones being reached by the lower limbs GB12265.2-2000
工业金属管道工程施工及验收规范 GB50235-2010	Code for construction and acceptance of Industrial Metallic Pipe Engineering GB50235-2010
电气装置安装工程电缆线路施工及验收规范 GB50168-2009	Code for construction and acceptance of cable system electric equipment installation engineering GB50168-2009
电气装置安装工程接地装置施工及验收规范 GB50169-2009	Code for construction and acceptance of earthing system electric equipment installation engineering GB50169-2009
机械电气安全 机械电气设备第一部分: 通用技术 EN60204-1-2009	Electrical safety of machine Mechanical Electrical Equipment Part 1:General Technology EN60204-1-2009
药品 GMP 实施指南 (2011 出版)	GMP Implementation Guide (2011)
WHO 数据与记录管理规范指南 2015.09	WHO Guidance on Good Data and Record Management Practices 2015.09
美国机械工程协会生物工艺设备指南, 2019 年 6 月 10 日	ASME BPE Bioprocessing Equipment, June 10, 2019

适用标准和法规	Applicable Standards And Regulations
ISPE 良好的自动化生产实践指南_遵从 GxP 计算机化系统监管的风险管理方法, 第 5 版	ISPE GAMP5 A Risk-Based Approach to Compliant GxP Computerized System, 5th Edition
ISPE 指南 卷 5: 调试与确认 (第一版)	ISPE Volume 5: Commissioning and Qualification (First Edition)


6 系统描述 SYSTEM DESCRIPTION

The scope of supply will cover four sets of stainless steel made purification modules, including 1 set of harvest tank, 3 sets of HPLC support functions, one set of depth filtration skid, 5 sets of purification tank (3×4000L, 2×3000L), UF support functions (three production lines, one set for each line; the fourth line, two sets), one set of 1500L gel tank(separate URS), one set of NF skid and bulk filling functions with their associated devices, pipeworks and accessories for two production lines. Each purification module, which can be further divided into pre-viral purification module and post-viral purification module, serves one dedicated production line. The four sets of purification modules are independent and physical segregated.

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The basic setup for one production line with main interfaces has been shown below.

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Below are some general practice for test data recording:

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Ensure that each document is “traceable” with a title, system reference (for example:Test No.).

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Each comments, entry, etc., on a document must be clear, readable, dated and signed. Do not use pencils to record data, but should use black or blue pen to record data.

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To verify the correct or incorrect information in a specific test, highlighters can be used. Test sheets with highlighted items shall include a description of the used colors and what they represent. The test sheets shall also follow the guidelines above regarding appendices and initial and date.

1) 绿色荧光– 相符合的信息/已安装 Green highlighter – consistent information/installed

2) 红色荧光– 不相符的信息/未安装 Red highlighter – disconsistent information/uninstalled

3) 黄色荧光–非关键的信息 Yellow highlighter – non critical information

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For each test executed, it is necessary to:

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Report the result, when it is applicable, in numerical form (it is not sufficient pass/fail if the result is a value)

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Note the procedure used for the test if it is not described in the body format used to execute the test

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Correct each verified document (wiring diagram, P&ID, layout, etc.) that is not corresponding to actual system, showing the as-built situation, these documents will be the base for future updating

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Attach each print-out/attachment produced during the test execution, signing and dating it, writing a clear Appendix No. and Test No. at which test the print-out/attachment refer to, if is present more than on attachment for the same test, Number progressively the attachment (page X of Y), signing and dating each attachment

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Sign and date each test (signatures of those who execute and verify the tests must be present)

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If some certain test is not executed, indicate the reason, and give, in the deviation form, a possible date for the execution

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The date format is yyyy-mm-dd, the format of signature must be consistent with the signature log .

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Do not use covering liquids or materials. Every writing mistake must be corrected as follows:

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Space not used must be barred and sign, or filled out with “NA”, which means “Not Applicable”.

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Deviation processing shall be carried out in accordance with Innovent’s internal procedures and the document “Deviation Management Procedure of Qualification”(SMP00209).

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Identify the information and qualification of all people that involved in protocol execution.

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List and identify all personnel involved in the execution of the present protocol (Name. Signature and Department/Company).

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Verify all personnel involved in the DQ execution have suitable qualification.

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All personnel involved in the execution of the present protocol should be appropriately trained DQ protocol in order to execute the protocol.

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All the operators are recorded by Name, Signature and Department/Company.

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All personnel involved in the DQ execution have suitable majority qualification.

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All the personnel who will execute this DQ protocol have been trained the protocol.

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可接受的标准 Acceptance Criteria			是 Yes /否 No
所有的参与人员已确认姓名、签名和部门/公司。 All the participants are identified by Name, Signature and Department/Company.			
参与确认的人员有相关的专业资质。 All personnel involved in the DQ execution have suitable majority qualification.			
所有执行本方案人员已接受培训。 All the personnel who will execute this DQ protocol have been trained the protocol.			
备注 Comments:			
偏差编号 Deviation No.			
执行人 Executed by		日期 Date	
确认人 Verified by		日期 Date	

9.2 设计需求文件的确认 DESIGN REQUIREMENT DOCUMENTS VERIFICATION

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Verify the supplied design requirement documents are available.

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Record the requirement documents generated from communications in project are available, include URS, meeting summary, mail and so on, then record document title, version and date of issued.

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All required documents must be available.

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9.3 设计文件的确认 DESIGN DOCUMENTS VERIFICATION

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Verify the design documents are available and standard.

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Record the design documents and drawings required by the DQ execution in the table below, and verify the document name, document number, version and the relevant approval status for each document.

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All the design documents and drawings required have been issued and available, the documents name, number and version have been recorded.

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Verify that design documents and drawings are compatible with URS.

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Input the each requirement of URS one by one in the table below.

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Identify the classification of each requirement which can be classified as Q, C, I. The criteria for the three categories are as follows:

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Category Q(Qualification): Critical requirements related to GMP, which will impact equipment performance, product quality, data integrity and patient safety. Need to respond that within design documents, and need to confirm during commissioning and qualification phase.

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Category C(Commissioning): Requirements for HSE and Non-GMP related requirements, which will not impact equipment performance, product quality, data integrity and patient safety. Need to confirm during commissioning phase.

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Category I(Instruction): other requirements for specification and information that won't involve design aspects, no need to respond. Such as project description, project delivery and warranty.

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Utilize the relevant specific information from design documents and drawings to respond every "C" requirements identified. The respond to the "N" requirements is to provide the reference only, and It doesn't need to respond the identified "I" requirements during DQ phase.

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The requirements in URS have been responded to by design documents, and the response is acceptable.

9.4.4 URS 符合性的确认检查表 URS COMPLIANCE VERIFICATION CHECK TABLE

M2 填料匀浆罐用户需求标准 M2 URS about Homogenize Tank System

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URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
					到 4000L 收获罐。	
URS4	<p>The temperature of the filtrated product will be cooled down and maintained at 2 ~ 8 °C in the 4000L harvest tank by chilled water in jacket side, for sample loading to the next step chromatography process.</p> <p>4000L 收获罐夹套会通入冷冻水使过滤后的产品降温并维持在 2-8°C，等待下步层析工艺。</p>	Q	<p>管道仪表流程图</p> <p>P&ID</p>	E01-203201-1-121	<p>The temperature of the filtrated product will be cooled down and maintained at 2 ~ 8 °C in the 4000L harvest tank by chilled water in jacket side, for sample loading to the next step chromatography process.</p> <p>4000L 收获罐夹套会通入冷冻水使过滤后的产品降温并维持在 2-8°C，等待下步层析工艺。</p>	
URS5	<p>The purification process will be carried out by 3×HPLC skids(out of scope), depth filtration skid, UF skids(out of scope) (three lines with one set of UF each line, the fourth line with two sets of UF) and NF skid to refine the target product from the supernatant solution of harvest process. Intermediates of purification process will be collected, held and adjusted in 5×purification tanks and disposable bags(out of scope).</p> <p>纯化过程从收获上清液开始，经过由 3 个层析模块（范围外）、深层过滤模块，UF 模块（范围外）(3 条线，每条线有 1 套 UF，第 4 条线有 2 套 UF)和 1 个纳滤模块达到纯化目标产品的目的。纯化过程中间体将由 5 个中间罐和一次性塑料袋（范围外）进行收集、存储和调配。</p>	Q	<p>管道仪表流程图</p> <p>P&ID</p>	E01-203201-1-121~133	<p>The purification process will be carried out by 3×HPLC skids(out of scope), depth filtration skid, UF skids(out of scope) (three lines with one set of UF each line, the fourth line with two sets of UF) and NF skid to refine the target product from the supernatant solution of harvest process. Intermediates of purification process will be collected, held and adjusted in 5×purification tanks and disposable bags(out of scope).</p>	

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
URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
	<p>line, finally processed by UF skid #2, and the concentrated product from UF operation will be transferred into disposable bag (out of scope) for bulk filling operation. For other three line, the product after NF and collected into 3000L purification tank without UF operation will be transferred into disposable bag (out of scope) for bulk filling operation.</p> <p>对病毒后的工艺，经过纳滤工艺后的纯化中间体将收集到 3000L 纯化罐。第四条生产线，产品进行超滤#2 的最终工艺操作，超滤操作的浓缩产物将被转移到一次性灌装袋（范围外）进行灌装操作。其他三条生产线，经过纳滤工艺后的纯化中间体将收集到 3000L 纯化罐，将直接被转移到一次性灌装袋（范围外）进行灌装操作。</p>				<p>器，也可以收集到 4000L 纯化罐，收集到混合器的中间体可以转移至 4000L 纯化罐，可以连接到深层过滤系统，也可以连接到 HPLC#2，收集到的中间体将在 4000L 纯化罐内进行低 pH 值的灭活。</p> <p>Depth filtration will be applied after low pH inactivation operation to remove impurities derived from process. Prior to depth filtration, intermediates in 4000L purification tank (TA03) should be processed for dilution, and the depth filtered product will be collected in 4000L purification tank(TA04), It can also be collected in a mixer.</p> <p>低 pH 值灭活后，经由深层过滤工艺来去除过程中产生的杂质。深层过滤之前，中间体在 4000L 纯化罐（TA03）进行稀释，过滤产品将在 4000L 纯化罐（TA04）进行收集，也可以收集在混合器里。</p> <p>After that, intermediates will</p>	


URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
					<p>be processed by HPLC #2 with elution collected in 4000L purification tank(TA05) ,It can also be collected in a mixer. Then you can chromatography HPLC#3 or UF/DF1 or NF.</p> <p>此后，中间体将进行层析 HPLC#2，洗脱液收集在 4000L 纯化罐（TA05），也可以收集在混合器里。随后可以进行层析 HPLC#3 或者 UF/DF1 或者 NF。</p> <p>Following adjusting of intermediates, HPLC #3 will be operated for polishing purification with elution collected in 3000L purification tank(TA06) ,It can also be collected in a mixer.</p> <p>随后进行中间体调整，并进行层析 HPLC#3 洗脱液收集在 3000L 纯化罐（TA06）也可以收集在混合器里。</p> <p>Then, intermediates will be further processed by UF skid #1 and NF skid with the processed product</p>	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
					<p>transferred to post-viral room, which is physically segregated with previous process detailed.</p> <p>然后，纯化中间体将进行超滤和纳滤工艺，产品转移到除病毒后房间进一步处理，除病毒前和除病毒后房间有明确的物理隔离。</p> <p>For post-viral process, intermediates from NF skid will be collected into 3000L purification tank. For the fourth line, finally processed by UF skid #2, and the concentrated product from UF operation will be transferred into disposable bag (out of scope) for bulk filling operation. For other three line, the product after NF and collected into 3000L purification tank without UF operation will be transferred into disposable bag (out of scope) for bulk filling operation.</p> <p>对病毒后的工艺，经过纳滤工艺后的纯化中间体将收集到 3000L 纯化罐。第四</p>	

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URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS16	 <p>The tank functions and piping requirements shall conform to the design shown above 管罐功能和管路需求应符合上图设计</p>	Q	管道仪表流程图 P&ID	E01-203201-1-127~133	<p>The tank function and pipeline design conform to the PFD. As shown in the P&ID.</p> <p>管罐功能和管路设计符合 PFD。详见 P&ID。</p>	
URS17	<p>The selections of the agitators for all tanks of purification module shall ensure the agitators can be ease of installation, maintenance, cleaning and are with good process performance. 所有纯化罐搅拌器的选择应保证搅拌器可以方便安装、维修、清洁，并具有良好的工艺性能。</p>	C	容器图纸 Vessel Drawing	WC02-26314 WR02-26315 WR02-26316 WR02-26328	<p>All purifier agitators can be easily installed, maintained, cleaned and have good process performance.</p> <p>所有纯化罐搅拌器可以方便安装、维修、清洁，并具有良好的工艺性能。</p>	
URS18	<p>Agitator for all tanks shall be variable speed, bottom-mounted sanitary magnetic agitator. The system shall consist of several elements – motor, VFD, gearbox, coupling, seal assembly, bearing assembly, shaft, and impeller(s). 所有罐的搅拌器应该是速度可调，底部安装的卫生磁力搅拌器。该搅拌器包括几个部件—电机、变频器、变速箱，耦合，密封组件、轴承组件、轴和叶轮。</p>	Q	磁力搅拌数据表 Magnetic Agitator Datasheet	PP20-0098-T-000379	<p>Agitators have variable speed, sanitary. The system shall consist of several elements – motor, VFD, gearbox, coupling, seal assembly, bearing assembly, shaft, and impeller(s). 磁力搅拌为卫生级设计，</p>	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
					可计速，磁力搅拌由电机、变频器、变速箱，耦合，密封组件、轴承组件、轴和叶轮组成。	
URS19	Agitator design of purification tank should ensure good mixing performance upon both the minimal working volume and the maximum working volume. 纯化罐搅拌器的设计应保证在最小的工作容积和最大工作容积都具有良好的混合性能。	Q	N/A	N/A	The mixing performance will be tested at FAT and SAT 将在 FAT,SAT 时测试搅拌的混合性能	
URS20	No split hubs will be allowed. Impeller components shall be machined or forged, not cast. One-piece welded construction of impeller is preferred whenever practical. Underside cavities in the hub, mechanical connections and crevices in the assembled wetted parts shall be minimized in order to facilitate cleaning-in-place and steam sterilization. Hubs shall be welded to the shaft unless noted otherwise. Welds shall be ground smooth to surface finish listed on data sheet. 非拆卸的轮轴是允许的。叶轮部件应是机械加工或锻造，不得铸造。实际应用时，叶轮首选整体焊接结构。为了便于就地清洗和蒸汽灭菌，应尽量减少轮毂中的下侧空腔、组装湿部的机械连接和缝隙，除非另有说明，应将毂焊接到轴上。焊缝表面光洁度应满足于数据表上列出的要求。	Q	磁力搅拌数据表 Magnetic Agitator Datasheet	PP20-0098-T-000379	The stirring shaft and the welding plate are not joined together. And for mechanical forging processing, CIP and SIP conveniently 搅拌轴和焊接盘不拼接。且为机械锻造加工，方便在线清洁和蒸汽灭菌。	
URS21	Supplier shall test agitator shaft for balance. Any vibration observed in factory testing and site testing at the owner's location shall be within industry standards. The supplier's shop drawings shall state first critical speed of each agitator. Mixing system shall be designed to avoid critical speed within normal operating range of the agitator. Any critical speed that falls within the operational speed range of the agitator shall be reported to the owner. Each agitator shall come up to full speed	C	N/A	N/A	Not Adaptable for Magnetic Agitator	

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URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS40	<p>To the maximum extent, the product wetting surface in purification system, including tanks with related transfer piping, shall be subjected to automated CIP and SIP procedure before or after each purification batch, or when necessary.</p> <p>尽可能地将纯化系统中产品润湿表面，包括带有相关输送管道的储罐，都必须在每次纯化生产前或之后进行自动 CIP 和 SIP 程序，或在必要时进行。</p>	Q	<p>功能说明</p> <p>FS</p>	<p>PP20-0098-T-004218~4222,</p> <p>PP20-0098-T-004206~4212,</p> <p>PP20-0098-T-009437</p>	<p>The product wetting surface in purification system, including tanks with related transfer piping, can be subjected to automated CIP and SIP procedure before or after each purification batch, or when necessary.</p> <p>纯化系统中产品润湿表面，包括带有相关输送管道的储罐，可以在每次纯化生产前或之后进行自动 CIP 和 SIP 程序，或在必要时进行。</p>	
URS41	<p>CIP solution will be supplied from CIP skid (not in purification module battery limit), and be pumped back to it or sent to drain after cleaning. One dedicated CIP skid will be provided for the pre-viral purification module in each production line. Another two CIP skids, which will be shared by buffer preparation and holding module (out of scope), will be provided for the post-viral purification module in each production line.</p> <p>CIP 溶液将从 CIP 模块（非纯化模块范围）供应，并在清洗后将其泵回或送至排水管。每个生产线的除病毒前纯化模块将使用一个专用 CIP 模块。缓冲液制备、储备模块和除病毒后纯化模块将共享 2 个 CIP 模块。</p>	Q	<p>管道仪表流程图</p> <p>P&ID</p>	<p>E01-203201-1-121~133</p>	<p>CIP solution will be supplied from CIP skid (not in purification module battery limit), and be pumped back to it or sent to drain after cleaning. One dedicated CIP skid will be provided for the pre-viral purification module in each production line. Another two CIP skids, which will be shared by buffer preparation and holding module (out of scope), will be provided for the post-viral purification module in each production line.</p> <p>CIP 溶液将从 CIP 模块（非</p>	

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URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
					于清洗。CIP 回水泵有低点排放，能完全排尽。	
URS54	<p>Centrifugal pump rated capacity shall not exceed capacity at the best efficiency point. Pump design shall permit installation of a larger impeller that would increase rated head by as much as 10%. Fabricate casings to allow removal and replacement of impellers (centrifugal).</p> <p>离心泵额定功率不应超过最佳效率点的功率。泵的设计应允许安装较大的叶轮，使额定扬程增加 10%。泵头保护套可拆卸，能更换桨叶。</p>	C	泵数据表 Pump Datasheet	PP20-0098-T-000082~83	<p>Centrifugal pump rated capacity does not exceed capacity at the best efficiency point. Pump design permit installation of a larger impeller that would increase rated head by as much as 10%. Fabricate casings to allow removal and replacement of impellers (centrifugal).</p> <p>森松选用离心泵额定功率不超过最佳效率点的功率。泵的设计允许安装较大的叶轮，使额定扬程增加 10%。泵头保护套可拆卸，能更换桨叶。</p>	
URS55	<p>Bearings shall be open, flush-through type (to prevent over-lubrication). Lubrication shall be with food-grade grease or oil. Pumps shall be furnished with o-ring seals. 轴承应打开，冲洗类型（防止过度润滑）。润滑应使用食品级润滑脂或机油。泵应该配置 O 形环密封。</p>	Q	泵数据表 Pump Datasheet	PP20-0098-T-000082~84	<p>Bearings can open, flush-through type. No need of lubrication oil product and water self-lubricate. Pumps is furnished with o-ring seals.</p> <p>轴承可打开冲洗类型。不需要润滑油产品和水自润滑。泵配置 O 形环密封。</p>	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS56	Flexible hoses shall be suitable for Compressed Air, Product, WFI, and caustic and acid CIP based applications. Hose shall be suitable for full vacuum operation including pump suction applications. Hose shall be autoclavable, constructed of stainless steel braid reinforced extruded platinum-cured silicone meeting USP XXI Class VI requirements. 软管适用于压缩空气、产品、注射用水以及 CIP 使用的碱和酸，SIP 应用。软管应适用于全真空操作，包括泵进口的应用。软管应高压灭菌，不锈钢编织增强铂金硫化硅胶的结构，满足 USP XXI VI 级要求。	Q	Flexible hose datasheet	PP20-0098-T-000911,000913,000914,000918~000920000925~000926	Meet the requirements. 满足。	
URS57	Outer layer shall have a rubber protective coating. Hose shall be of smooth bore construction and shall be manufactured for sanitary applications. Tubing fittings shall be factory installed, integral with the hose, and shall be of the sanitary clamp type. Length of hose shall be as required for application. 外层应有橡胶保护层。软管应具有光滑的孔结构，可应用于卫生场合。管道配件应在工厂安装，与软管整体连接，并应采用卫生卡接。软管长度应按要求使用。	Q	Flexible hose datasheet	PP20-0098-T-000911,000913,000914,000918~000920000925~000926	Meet the requirements. 满足。	
URS58	For better sanitary performance, the supplier shall design the skid avoid of flexible hoses as much as possible. 为了达到更好的卫生性能，供应商在模块设计时应尽可能避免使用软管。	C	管道仪表流程图 P&ID	E01-203201-1-121~133	The least flexible hoses are designed in PID.PID 中尽量少的设置软管。	
URS59	All the tanks shall be furnished with appropriate tank bottom valve. The valve shall be self-drainable. The valve body can be automatically CIP'ed and SIP'ed through dedicated port without impact content in tanks. 所有的罐都配有适当的罐底阀，且阀门为自排净设计。阀体能够自动的进行 CIP/SIP，而不影响罐内。	Q	管道仪表流程图 P&ID 罐底阀数据表	E01-203201-1-121,123,125,127,130,132 PP20-0098-T-000438 PP20-0098-T-000443	All tanks are equipped with plunger - type tank bottom valves The valve shall be self-drainable. The valve body can be automatically CIP'ed and SIP'ed through dedicated port without impact content in tanks. 所有的罐配备柱塞式罐底阀，阀门为自排净设计。阀体能够自动的进行	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
					CIP/SIP, 而不影响罐内。	
URS60	<p>For all tanks, disposable sampling valve which is pre-sterilized shall be provided. Sampling bag which matches sampling valve shall be provided by client. The operation of the sampling valve shall have no impact on the product inside the tank.</p> <p>对于所有的罐, 使用预消毒的一次性取样阀, 业主自备配套的取样袋。取样阀的操作应对罐内产品无影响。</p>	Q	取样阀数据表 Sample Valve Datasheet	PP20-0098-T-000440 PP20-0098-T-000831	<p>All vessels are sampled with pre-sterilized disposable sampling bags provided by the owner. The operation of the sampling valve has no effect on the product in the tank.</p> <p>对于所有的罐, 使用预消毒的一次性取样袋取样, 业主自备配套的取样袋。取样阀的操作对罐内产品无影响。</p>	
URS61	<p>Using pure steam control valve, the valve opening can be automatically adjusted according to the pressure setting value to meet different steam consumption requirements. Valve shall have appropriate design pressure and temperature.</p> <p>选用纯蒸汽调节阀, 可根据压力设定值, 自动调节阀门开度, 满足不同蒸汽用量需求。阀门应该合理设计压力和温度。</p>	C	调节阀数据表 Control Valve Datasheet	PP20-0098-T-000850	<p>Using pure steam control valve, the valve opening can be automatically adjusted according to the pressure setting value to meet different steam consumption requirements. Valve have appropriate design pressure and temperature.</p> <p>选用纯蒸汽调节阀, 可根据压力设定值, 自动调节阀门开度, 满足不同蒸汽用量需求。阀门压力和温度合理。</p>	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS62	Provide self-adjusting balanced thermostatic type steam trap designed for clean steam sanitary operation. 为洁净蒸汽灭菌操作提供自平衡式热静力蒸汽疏水阀。	C	疏水阀数据表 Steam Trap Datasheet	PP20-0098-T-000853~854	Provide self-adjusting balanced thermostatic type steam trap designed for clean steam sanitary operation. 森松为洁净蒸汽灭菌操作提供自平衡式热静力蒸汽疏水阀。	
URS63	Process systems shall be designed so that the risk of contamination/cross-contamination is minimized between media that shall not get in contact with each other. 工艺系统的设计应使介质之间的污染/交叉污染的风险最小化。	Q	管道仪表流程图 P&ID	E01-203201-1-121~136	Process systems are designed so that the risk of contamination/cross-contamination is minimized between media that shall not get in contact with each other. 工艺系统的设计满足介质之间的污染/交叉污染的风险最小化。	
URS64	Prevention against contamination/cross-contamination through leaking valves shall be properly established between CIP systems and other media and always between water systems and other media. 在 CIP 系统和其他介质之间，水系统和其他介质之间适当设置防止污染/交叉污染的阀门。	Q	管道仪表流程图 P&ID	E01-203201-1-121~136	Prevention against contamination/cross-contamination through leaking valves are properly established between CIP systems and other media and always between water systems and other media. 在 CIP 系统和其他介质之间，水系统和其他介质之间适当设置防止污染/交叉	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
					污染的阀门。	
URS65	Drainage towards drains shall be secured against reverse suction and contamination with air breaks. Alternatively, a suitable sanitary mechanical device may be used, if the drain connection needs to be closed. 排水管道应防止反向吸入和空气污染。如果需要关闭排水连接，则可以使用合适的卫生型机械装置。	Q	管道仪表流程图 P&ID	E01-203201-1-121~136	Air break drains are installed at the customer docking position and check valve are installed at our necessary position. 在信达对接位置设置防倒灌地漏，在我方必要位置设置止回阀。	.
URS66	All the instruments on SF1/SF4 shall be of sanitary design. Instruments installed within the sterile envelope or boundary shall be designed for SIP. Consideration shall be made in the design for instrument removal for calibration. 管道等级为 SF1/SF4 的仪表均应卫生设计。安装在无菌区域或边界内的仪器应设计为 SIP。设计时应考虑仪表的拆卸和校准。	Q	压力变送器、流量计、压力表数据表 PRESSURE TRANSMITTER 、 MASS FLOWMETER 、 PRESSURE GAUGE Datasheet	PP20-0098-T-000461 PP20-0098-T-000865 PP20-0098-T-000465 PP20-0098-T-000162 PP20-0098-T-000845 PP20-0098-T-000961	All the instruments on SF1/SF4 are of sanitary design. Instruments installed within the sterile envelope or boundary are designed for SIP. Consideration are made in the design for instrument removal for calibration. 管道等级为 SF1/SF4 的仪表应卫生设计。安装在无菌区域或边界内的仪器设计为 SIP。设计时考虑了仪表的拆卸和校准。	
URS67	Instruments installed within the sterile envelope or boundary shall be designed for CIP or removed for COP. In the case of COP, blind caps or plugs shall be provided to maintain the integrity of the system. 安装在无菌区域或边界内的仪器应设计为 CIP 或 COP。在离线清洗的情况下，应提供盲板或堵头以保	Q	管道仪表流程图 P&ID	E01-203201-1-121~136	Instruments installed within the sterile envelope or boundary design for CIP. In the case of COP, blind caps or plugs are provided to maintain the integrity of the	

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URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
					提供具有传感器管理系统的智能 pH 传感器, 可以测量 pH 传感器的寿命、CIP / SIP 循环、校准历史等。	
URS71	<p>Level / Weight Instrumentation System 液位/称重测量系统</p> <p>Load cell or scale or level sensors (indicated in PIDs) shall be used for weighing of all tanks or bottles or bags. The installation of load cell and fabrication of supporting platform (if applicable) shall comply with the installation instructions from the load cell supplier.</p> <p>称重模块或秤或液位传感器 (显示在 PID) 用于所有罐或瓶或袋称重。称重模块的安装和支撑平台的制作 (如适用) 应符合称重模块供应商的安装说明。</p> <p>The weighing accuracy of all scales/balances shall be less than $\pm 20g$.</p> <p>所有的秤/天平称量精度应小于$\pm 20g$。</p> <p>The design of balance / scale shall be ease of cleaning and sanitization, shall be suitable for operation requirements in clean room.</p> <p>秤/天平的设计应易于清洗和消毒, 应适用于在洁净室的操作要求。</p>	Q	称重数据表	PP20-0098-T-000448 PP20-0098-T-001026	<p>Load cell or scale or level sensors (indicated in PIDs) shall be used for weighing of all tanks or bottles or bags. The installation of load cell and fabrication of supporting platform (if applicable) can comply with the installation instructions from the load cell supplier.</p> <p>称重模块或秤或液位传感器 (显示在 PID) 用于所有罐或瓶或袋称重。称重模块的安装和支撑平台的制作 (如适用) 符合称重模块供应商的安装说明。</p> <p>The weighing accuracy of all scales/balances are $\pm 10g$.</p> <p>所有的秤/天平称量精度为$\pm 10g$。</p> <p>The design of balance / scale can be ease of cleaning and sanitization, shall be suitable for operation requirements in clean room.</p> <p>秤/天平的设计易于清洗和消毒, 应适用于在洁净室的操作要求。</p>	

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URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS93	<p>The surface finish requirement shall be compliant with general piping specifications and vessel datasheet. 表面光洁度要求应符合管道等级表和罐体配置清单（见 URS 附件）。</p> <p>Internal surfaces in contact directly with product and WFI supply line shall be SF4 内表面直接接触产品和注射用水的管线为 SF4 Clean steam post regulator shall be SF1 at least. 清洁蒸汽调压阀后应至少为 SF1。 CIP supply and return lines shall be SF1 CIP 供应和回流管道应为 SF1 Drain line, kill drain line, and clean condensate post block valve shall be SS304 排放管线，杀菌排放管线，和清洁冷凝水切断阀后应为 SS304 不锈钢 Process gases and vent lines post sterile filters shall be SF1 at least. 工艺空气和无菌过滤器后排气线应至少为 SF1。</p>	Q	Piping Class Index 管道等级索引表	PP20-0098-T-000012	<p>The surface finish requirement will be compliant with general piping specifications and vessel datasheet. 表面光洁度要求符合管道等级表和罐体配置清单（见 URS 附件）。</p> <p>Internal surfaces in contact directly with product and WFI supply line is SF4 内表面直接接触产品和注射用水的管线为 SF4 Clean steam post regulator shall be SF1 . 清洁蒸汽调压阀后为 SF1。 CIP supply and return lines is SF1 CIP 供应和回流管道为 SF1 Drain line, kill drain line, and clean condensate post block valve is SS304 排放管线，杀菌排放管线，和清洁冷凝水切断阀后为 SS304 不锈钢 Process gases and vent lines post sterile filters is SF1 . 工艺空气和无菌过滤器后为 SF1。</p>	
URS94	<p>Other lines and equipment which connecting with black utilities and without direct contact with product will be figured out by supplier to comply to current industry standard and operation requirements. 其他与公用设施连接且与产品无直接接触的线路和设备将由供应商设置，以符合现行的行业标准和操作要求。</p>	C	Piping Class Index 管道等级索引表	PP20-0098-T-000012	<p>Other lines and equipment which connecting with black utilities and without direct contact with product is figured out by SMP to comply to current industry standard and operation requirements.</p>	

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URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS97	Supplier shall propose method to be used to assure weld quality. Client to review controls before manufacturing commences. Minimum 20% of all product/WFI contact welds shall be tested/inspected and certified. 供货商应提出措施来保证焊接质量。在开始制造前，客户将确立审查控制。20%以上所有产品/注射用水接触的焊缝应测试/检验。	Q	Piping Class Index 管道等级索引表	PP20-0098-T-000012	20% orbital welds require borescope inspection	
URS98	All weld beam shall be correctly tagged and tag number shall be permanently marked on pipe in correspondence with that in documents. 所有焊缝应正确标记和标签号码应永久标记在管道并与相关文件对应。	Q	Piping Design Specification 管道设计说明	PP20-0098-T-010289	All welds will be correctly tagged and tag number will be permanently marked on pipe in correspondence with that in documents. 所有焊缝将正确标识，同时标识号码永久地标记在管路上，并与在文件存档的内容一致。	
URS99	The design shall aim at including as few deadlegs as possible. The identification of dead leg ("min." in P&D) is "3D rule". 设计应尽可能减少死角。死角的标识 ("min") 是 "3D" 规则。	Q	Piping Design Specification 管道设计说明	PP20-0098-T-010289	The means "3D rule" in P&ID by "min", Piping design followed this rule.	
URS100	Short pipe sections shall preferably be designed with a 1% slope and long pipe sections shall be designed with a 0.5% slope. Slopes below 0.5% can only be accepted in exceptional cases. 短管段最好设计成 1% 坡度，长管段应采用 0.5% 坡度设计。低于 0.5% 的斜坡只能在特殊情况下接受。	Q	Piping Class Index 管道等级索引表	PP20-0098-T-000012	Short pipe sections is designed with a 1% slope and long pipe sections is designed with a 0.5% slope. Slopes below 0.5% can only be accepted in exceptional cases. 短管段采用 1% 坡度，长管段采用 0.5% 坡度设计。	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS101	Equipment shall be designed to be easy to clean. 设备应设计成易于清洗。	I	容 器 图 纸 Vessel Drawing	WC02-26314 WR02-26315 WR02-26316 WR02-26328	The equipment is designed according to sanitary grade and easy to clean 设备按照卫生级设计，易于清洗	
URS102	Sharp corners, edges or nuts that tear cleaning tissues and prevent easy wiping are not allowed 不允许出现会撕破擦拭工具的尖角，锋利的边缘以及不易擦拭的地方。	C	Piping Design Specification 管道设计说明	PP20-0098-T-010289	Module sharp corners, edges will be polished or chamfered, which will make the cleaning easily. 模块尖角、边缘都进行打磨或倒角处理，易于清洁。	
URS103	The design of all equipment that will contain pharmaceutical materials shall ensure that all water will be self-drain, to achieve efficient drying. 所有含有制药原料的设备的设计应确保所有的水都能自排净，以实现高效干燥。	Q	管道仪表流程图 P&ID	E01-203201-1-121~136	Pipe are desinged with slope to drain. 管道设计带有坡度可使系统排尽。	
URS104	Parts that hide or trap spilled product are to be avoided. 避免出现会藏留产品的部位。	Q	管道仪表流程图 P&ID	E01-203201-1-121~136	The design is aiming at min deadleg. 设计已尽可能减少死角。	
URS105	All the materials used in cleanroom shall be compatible with sanitizers, such as IPA, sporicidal agent. 所有用于洁净室材料应能耐受消毒剂，如 IPA，杀菌剂。	C	Piping Class Index 管道等级索引表 Piping Design Specification 管道设计说明	PP20-0098-T-000012 PP20-0098-T-010289	All piping compnents and valves materials used in cleanroom are compatible with sanitizers, such as IPA, sporicidal agent. 所有用于洁净室管道阀门材料能耐受消毒剂，如 IPA，杀菌剂。	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS106	<p>All bearings and mechanical seals shall be designed to prevent particulate matter from being shed into the process or direct utilities. Bearings shall be stainless steel and sealed-for-life</p> <p>所有轴承和机械密封的设计应能防止颗粒物质洒到或直接接触轴承内部。轴承应为不锈钢和永久性密封。</p>	Q	<p>Magnetic Agitator Datasheet</p> <p>磁力搅拌数据表</p> <p>Pump Datasheet</p> <p>泵数据表</p>	<p>PP20-0098-T-000376</p> <p>PP20-0098-T-000381</p> <p>PP20-0098-T-000377</p> <p>PP20-0098-T-000082</p> <p>PP20-0098-T-000083</p> <p>PP20-0098-T-000084</p> <p>PP20-0098-T-000085</p> <p>PP20-0098-T-000062</p> <p>PP20-0098-T-000781</p>	<p>Bearings and mechanical seals of the magnetic stirring and pump are designed to prevent particulate matter from being shed into the process or direct utilities. Bearings shall be stainless steel and sealed-for-life</p> <p>磁力搅拌和泵轴承和机械密封能避免散落物质到工艺或公用工程中。轴承材质采用不锈钢并完全密封。</p>	
URS107	<p>All moving parts shall be lubricated to ensure maintenance free operation. The design of the lubrication equipment shall ensure that the lubricant cannot leak into the process, direct utilities or into the cGMP area.</p> <p>所有运动部件应润滑以保证免维护操作。润滑设备的设计应确保润滑油不会泄漏到工艺中，进入公用工程或 cGMP 的区域。</p>	C	<p>Magnetic Agitator Datasheet</p> <p>磁力搅拌数据表</p> <p>Pump Datasheet</p> <p>泵数据表</p>	<p>PP20-0098-T-000376</p> <p>PP20-0098-T-000381</p> <p>PP20-0098-T-000377</p> <p>PP20-0098-T-000062</p> <p>PP20-0098-T-</p>	<p>Magnetic stirring lubricants are designed to be maintenance-free and will not leak into the process and into utility or cGMP areas. Pump use product and water self-lubricate.</p> <p>磁力搅拌的润滑油为免维护的设计，泵是利用产品和水自润滑。润滑油不会</p>	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
				000781 PP20-0098-T-000446	泄漏到工艺中，进入公用工程或 cGMP 的区域	
URS108	Food grade lubricants that are easy to clean using WFI, shall be used for all locations where mechanical failure of lubricant seals would allow the lubricant to leak into areas where it may come into contact with the product. The supplier shall provide supporting documentation for each lubricant and use of lubricants will require prior approval by Owner. 易于用 WFI 进行清洁的食品级润滑油，将用于密封失效会让润滑油泄漏到可能接触到产品的区域的所有机械密封的润滑。供应商应提供润滑油支持文件，使用润滑的位置都需要事先由业主批准。	C	N/A	N/A	There is no lubricants that are potentially contacting products. 无润滑油会接触/潜在接触产品	
URS109	Pressure relief systems shall only be used where required by Chinese Standards where pressure sensitive or the pressure rating of downstream equipment requires protection. 泄压系统只能用于中国标准，压力敏感或下游设备的压力等级要求的防护要求。	C	管道仪表流程图 P&ID	E01-203201-1-121~136	Morimatsu designs rupture disc on vessl safety valves at the jacket outlet and the TCU industrial steam inlet to protect the container jacket and the TCU heat exchanger. 森松在夹套出口和 TCU 工业蒸汽入口设计安全阀，罐顶设置爆破片，用来保护容器夹套和 TCU 换热器。	
URS110	Pressure relief valve shall vent directly to atmosphere. 泄压阀应直接向大气排气。 Shall be suitable for operation over the design temperature and pressure ranges of the system. 应在系统设计的温度和压力范围内操作。	C	管道仪表流程图 P&ID 安全阀数据表 Pressure Safety	E01-203201-1-121~136 PP20-0098-T-001045~1046	Mori designed relief valve outlet piping to discharge directly to the surface. 森松设计安全阀出口管道直接排放到地面上。	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
			Valve Datasheet		Pressure safety valve and rupture disc can release pressure when system over pressure. 安全阀和爆破片可在系统超压时工作释放压力。	
URS111	Shall be installed as per the manufacturer's recommendations. 应按照制造商的建议安装。	C	安全阀数据表 Pressure Safety Valve Datasheet	PP20-0098-T-001045~1046	Provide valve manual. 提供阀门安装手册	
URS112	PR shall be used where required for compressed gases and steam utilities where protection to equipment and / or personnel is required. 在需要保护设备和/或人员的地方使用压缩气体和蒸汽公用设施的，压力调节阀应该被使用。	Q	管道仪表流程图 P&ID	E01-203201-1-121~136	The pressure-reducing valve is set on the compressed air pipeline, and the pressure control valve is set on the pure steam pipeline. 压缩空气管道上设置减压阀，纯蒸汽管路上设置压力控制阀。	
URS113	Sanitary types shall be 316 SS diaphragm balanced spring type with internal finishes as per specifications for sanitary equipment. 卫生类型应为 316 SS 膜片平衡弹簧类型，内部按照卫生设备的规格完成。	Q	Control Valve Datasheet 调节阀数据表	PP20-0098-T-000070 PP20-0098-T-000850 PP20-0098-T-004187 PP20-0098-T-004188	Sanitary control valves are made of SS316L diaphragm balanced spring type with Ra<=0.4um. 卫生型调节阀均为 SS316L 隔膜弹簧平衡型，内表面抛光为 Ra<=0.4um，满足无菌设备要求	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS114	Shall be installed as per the manufacturer's recommendations. 应按照制造商的建议安装。	C	Control Valve 调节阀数据表	PP20-0098-T-000070 PP20-0098-T-000850 PP20-0098-T-004187 PP20-0098-T-004188 PP20-0098-T-001042 PP20-0098-T-001044	Installation manual.is provided 提供安装手册	
URS115	Each instrument, component and valve, which is identified by a tag or equipment number on the P&ID or data sheets shall be supplied with a nameplate, engraved with the tag or equipment number and mechanically fixed to the equipment. The nameplate shall be of stainless steel. The information and format on nameplate shall be confirmed by owner's QA department before fabrication. 每一个仪表、部件和阀门，由标签或设备号上的 P&ID 或数据表标识，应提供一个铭牌，上面刻有标签或设备号，并机械地固定在设备上。铭牌应该是不锈钢的。铭牌上的信息和格式在生产前由业主 QA 部门确认。	Q	P&ID 管道仪表流程图	E01-203201-1-121-136	Each instrument, component and valve is identified by a tag number on the P&ID. The information and format on nameplate will be confirmed by owner's QA department during commissioning. P&ID 中体现了各部件的位号。铭牌的安装和正确性将在调试与确认活动中核查。	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS116	Generally the whole control system design will be in the scope of another dedicated automation system supplier while not in the scope of purification module fabrication supplier referred in this specification. However, the purification module fabrication supplier shall take responsibility on specification/selection, procurement, installation of all instruments which have been illustrated in owner's P&ID. 一般来说，整个控制系统设计将在另一个专用自动化系统供应商的范围内，而不在本规范中提到的纯化模块制造供应商的范围内。然而，纯化模块制造供应商应负责在 P&ID 范围内所有仪表的技术要求/选择，采购和安装所有仪器。	I	N/A	N/A	N/A	
URS117	The automation system supplier will take responsibility for wiring of all solenoids, instruments, motors and provide proper electrical supply to VFDs. 自动化系统供应商将负责所有的电磁阀，仪表，电机的接线工作并为变频器提供电源供应。	I	Charm IO 和 MCC 电气图、布局图 Electrical 、Layout Drawing of Charm IO and MCC	3287830-DWG- AC-SA-96NIS MCC-DS1-CIP-001	The automation system supplier can take responsibility for wiring of all solenoids, instruments, motors and provide proper electrical supply to VFDs. 自动化系统供应商将负责所有的电磁阀，仪表，电机的接线工作并为变频器提供电源供应。	
URS118	For automation interface requirements on components and instruments, refer to <i>Instrument General Specification</i> . 有关元件和仪器的自动化接口要求，请参阅仪表通用规范。	I	N/A	N/A	N/A	
URS119	Valve terminals and weighing system use Profibus DP. 阀岛和称重使用 Profibus DP 通讯。	C	称重数据表、控制架构图	PP20-0098-T-000448 PP20-0098-T-001026 3287830-PAS-DV1-DWG	Valve terminals and weighing system use Profibus DP. 阀岛和称重使用 Profibus DP 通讯。	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS120	Valve terminals is required that the supplier shall supply the filed panels with closure rating IP55. Panels are fixed up inside skid or nearby tanks. The longest distance from probe, valve, elements etc. shall not exceed 5 meters. 阀岛安装在箱子内，至少 IP55 等级。控制柜安装在模块上或罐体旁边。探头，阀门和部件之间的最长距离不超过 5 米。	C	阀岛箱布局图	PP20-0098-T-010241 PP20-0098-T-010242	Valve terminals is required that the supplier shall supply the filed panels with closure rating IP55. Panels are fixed up inside skid or nearby tanks. The longest distance from probe, valve, elements etc. shall not exceed 5 meters. 阀岛安装在箱子内，至少 IP55 等级。控制柜安装在模块上或罐体旁边。探头，阀门和部件之间的最长距离不超过 5 米。	
URS121	The pneumatic actuated valves solenoid banks are scope of modular system supplier. The bank can be either 8 units simplex type with Profibus connection. 气动阀阀岛是系统供应商的范围。阀岛是 8 一组，通过 Profibus 通讯。	C	阀岛 BOM 清单	PP20-0098-T-009858	The pneumatic actuated valves solenoid banks are scope of modular system supplier. The bank can be either 8 units simplex type with Profibus connection.气动阀阀岛是系统供应商的范围。阀岛是 8 一组，通过 Profibus 通讯。	
URS122	3Ph, 380V, 50Hz 3Ph, 380V, 50Hz 1Ph, 220V, 50Hz 1Ph, 220V, 50Hz Power supply shall be brought by owner to each main power cabinets only. All cables between control cabinets for machines shall be supplied by the supplier. 电源将由业主送到每个主电源柜。所有机器的控制柜之间的电缆应由供应商提供（应为自控供应商）。	C	MCC 柜图纸 The Drawing of MCC	MCC-DS1-PUR-001	3Ph, 380V, 50Hz 3Ph, 380V, 50Hz 1Ph, 220V, 50Hz 1Ph, 220V, 50Hz Power supply was brought by owner to each main power cabinets only. All cables between control cabinets for machines is supplied by the supplier. 电源由业主送到每个主电源柜。所有机器的控制柜之间的电缆由供应商提供（为自控供应商）。	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS123	Cable tray inside the skids shall be provided by supplier. The cables inside the scope of supplier shall be laid down by supplier. 模块内的电缆桥架应由供应商提供。供应商的范围内的电缆应由供应商敷设。	C	N/A	N/A	Cable tray inside the skids is provided by supplier. The cables inside the scope of supplier is laid down by supplier. 模块内的电缆桥架由森松提供。森松的范围内的电缆由森松敷设。	
URS124	All cables shall be LSZH flame retardant cable 所有电缆应无卤低烟阻燃电缆	C	电缆表 Cable List	PP20-0098-T-009646 PP20-0098-T-000826	All cables are LSZH flame retardant cable 所有电缆是无卤低烟阻燃电缆	
URS125	The use of connectors shall be maximizing with the advantage of reduced installation time by having a reliable connection on the equipment. 利用在设备上有可靠连接的连接器，其使用应最大程度地减少安装时间。	C	N/A	N/A	The use of connectors are maximizing with the advantage of reduced installation time by having a reliable connection on the equipment. 利用在设备上是有可靠连接的连接器，其使用最大程度地减少安装时间。	
URS126	All the cabinets, cables, wiring and conduit will be designed and installed by third party. The purification module supplier shall leave space and coordinate with third party during shop design and manufacturing. 所有的机柜、电缆、布线和管道都将由第三方设计和安装。纯化模块供应商应在车间设计和制造过程中留出空间并与第三方协调。	C	N/A	N/A	All the cabinets, cables, wiring and conduit is designed and installed by third party. The purification module supplier already leave space and coordinate with third party during shop design and manufacturing. 所有的机柜、电缆、布线和管道都将由第三方设计和安装。纯化模块供应商	

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URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
					Interconnecting terminals. 互连终端。	
URS130	<p>Live parts of electrical equipment (e.g. junction boxes, panel boxes, receptacle outlets, equipment cabinets, etc.) and/or junction points shall be provided with proper covers. Panels of cabinets shall be designed and built to meet the specified area classification. Unless otherwise specified all electrical cabinets and panels shall be IP54 or better.</p> <p>带电设备（如接线盒、面板盒、插座、设备柜等）和/或连接点的带电部分应有适当的覆盖物。柜面板应设计和建造，以满足指定区域的分类。除非另有说明，所有电气机柜及面板防护等级应满足 IP54 或更高。</p>	C	MCC 、 Charm IO Panel 柜图纸	MCC-DS1-PUR-001 3287830-DWG-AC-SA-96NIS	<p>Live parts of electrical equipment (e.g. junction boxes, panel boxes, receptacle outlets, equipment cabinets, etc.) and/or junction points is provided with proper covers. Panels of cabinets is designed and built to meet the specified area classification. Unless otherwise specified all electrical cabinets and panels is IP54 or better.</p> <p>带电设备（如接线盒、面板盒、插座、设备柜等）和/或连接点的带电部分有适当的覆盖物。柜面板应设计和建造，以满足指定区域的分类。除非另有说明，所有电气机柜及面板防护等级满足 IP54 或更高。</p>	
URS131	<p>Identify both ends of cables using slide on plastic markers with the same identification as shown on the wiring diagrams and terminals. Enclose cabling in slotted plastic trunking to 60% capacity of trunking only.</p> <p>接线两端的标识采用可滑动的塑料标签，标签名与接</p>	C	N/A	N/A	will be checked during FAT, SAT 将在 FAT,SAT 时核查	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
	线图和接线端子上的一致。导线槽的布线量是设计容量的 60%。					
URS132	Termination to external equipment must be via terminals. Connect all incoming and outgoing mA, mV and 24V signals to disconnect terminals type. Connect no more than two cores to one side of any terminal. Terminals shall be identified with the same identity as shown on the wiring diagrams. 通过接线端子可以中止与外部设备的通讯。连接所有进出的 mA, mV 和 24V 信号到断开型端子。任何端子一侧的连接不超过两根线。端子标识与接线图上一样。	C	MCC、Charm IO 柜 图纸	MCC-DS1-PUR-001 3287830-DWG-AC-SA-96NIS	Termination to external equipment must be via terminals. Connect all incoming and outgoing mA, mV and 24V signals to disconnect terminals type. Connect no more than two cores to one side of any terminal. 通过接线端子可以中止与外部设备的通讯。连接所有进出的 mA, mV 和 24V 信号到断开型端子。任何端子一侧的连接不超过两根线。	
URS133	The cable tray in skid shall be closed type and ease of cleaning. 模块内电缆桥架应封闭，便于清洁。	C	N/A	N/A	The cable tray in skid is closed type and ease of cleaning. 模块内电缆桥架封闭。	
URS134	Ensure all equipment items are fully earth bonded. 确保所有设备都是完全接地的。	C	N/A	N/A	All equipment items are fully earth bonded. 所有设备都是完全接地的。	
URS135	Provide an M6 earth stud on the inside and outside of metallic control panels and junction boxes for connection to the main earth system. Provide an earth bar inside control panel insulated from the panel, for signal earths and screens of external cables	C	N/A	N/A	will be checked during FAT, SAT 将在 FAT,SAT 时核查	

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URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS139	All elements, which are required for operational reasons, must be easily accessible from the production room and comply with cGMP. 所有由于操作原因所必需的元件，必须从生产房间容易进入，并遵守 cGMP。	C	Second floor Equipment Module Layout 2 楼设备模块布置图	PP20-0098-T-000021	All elements, which are required for operational reasons, are easily accessible from the production room and comply with cGMP. 所有由于操作原因所必需的元件，容易从生产房间进入，并遵守 cGMP。	
URS140	External surface temperature with insulation shall be lower than 45°C during all phases. Hot spots access shall be signaled and protected. 在所有阶段，外部表面绝缘温度应低于 45°C。热点接入应标识并加以保护。	C	N/A	N/A	will be checked during FAT, SAT 将在 FAT,SAT 时核查	
URS141	All equipment shall conform to an average of <75 dBA weighted scale and peak noise of 80 dBA when measured at 1 meter in any direction from the equipment when operating at all maximum loads and speeds. 所有设备在最大负荷和速度下运行时，在任何方向上 1 米范围外测量应满足平均小于 75 分贝，峰值 80 分贝噪声。	C	N/A	N/A	will be checked during FAT, SAT 将在 FAT,SAT 时核查	
URS142	P&ID change approval P&ID 变更和批准 The owner approved P&ID for construction will be provided. The P&ID shall be updated if changes occurred during the execution of the project. The final as-built P&ID will be sent to owner for approval prior FAT. 供应商负责 P&ID 图的详细设计，建筑 P&ID 图经业主批准后才可施工。在项目执行过程中，P&ID 图应及时更新。项目 FAT 前，最终版 P&ID 图应获得业主批准。	I	N/A	N/A	N/A	
URS143	Supplier shall design and deliver the process documents, eg Tag List, Datasheet, Calculation Sheets etc based on the provided P&IDs. 供应商应设计和提交工艺文件，如设备位号列表，设备技术文件，能耗计算书等列在交付清单上的文件。	I	N/A	N/A	N/A	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS144	<p>Datasheet verification 数据表确认</p> <p>Templates of ISO and ISA standardized datasheet are provided by owner. The supplier shall design datasheet for verifying all individual component and instrumentation shown on P&IDs.</p> <p>业主会提供 ISO 和 ISA 数据标准模板。收获系统供应商应提供设备数据列表用来查证 P&ID 上出现的每个组件和仪表。</p> <p>In a cost or time effective issue, Supplier may change to use different sub- supplier, or model of equipment. However, specification has to be approved by owner or its Engineering Consultant prior to use.</p> <p>考虑到设计费用和工期, 供应商也可以采纳二级供应商的资料和设备模型。但这些设备说明文件在使用前应经业主或第三方工程顾问的批准。</p> <p>Moreover, instrumentation datasheet shall be approved by owner appointed process automation system supplier to make sure the connectability of electrical signal and IA, as well as the effective control abilities on them. However, process fluidic and sanitary characteristics shall be verified by supplier.</p> <p>而且, 仪表数据列表应得到业主指定自控系统供应商的批准, 以确保仪表电气信号的连接和控制性能。然后工艺流体性能和卫生要求由收获系统供应商确认。</p>	I	N/A	N/A	N/A	
URS145	<p>Risk analysis 风险分析</p> <p>Risk analysis on the systems within scope of work (Failure mode and effects analysis, FMEA) 供货范围内的系统风险分析 (采用 FMEA: 失效模式和影响分析)</p>	I	N/A	N/A	N/A	
URS146	<p>Design Qualification (DQ) 设计确认</p> <p>Prior to the performance of fabrication, the supplier must complete the design qualification. The supplier shall develop a design qualification protocol for owner's review and approval before execution.</p> <p>在制造前必须完成 DQ。供应商应发布设计确认方案, 以供页数审核和批准。</p>	I	N/A	N/A	N/A	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS147	<p>GA release and approval 系统布置图</p> <p>Prior to the fabrication, the fabrication level GA (General Arrangement of system) in 3D model shall be designed by the supplier based on the P&IDs provided by owner. 在制造开始之前，供应商必须完成整个 3D 设计，并得到业主批准。在 3D 模型中英包含所有 P&ID 中管路和配件。</p> <p>The 3D model shall show actual equipment / components dimension, connection type and location etc for final fabrication with compliance to the dedicated equipment / component installation requirements. 供应商设计的带有 3D 示意，设备尺寸，管路连接类型和连接位置等设备系统布置文件和容器草图应得到业主的批准。</p> <p>3D GA shall be approved by owner, and the Vessel Drawing as well.</p> <p>3D 布置图和罐体设计图纸都应得到业主的批准</p>	I	N/A	N/A	N/A	
URS148	<p>Quality Plan 质量计划</p> <p>The supplier shall develop a Quality Plan for owner's review and approval. 供应商应提供一份产品质量计划给业主审阅和批准。</p>	I	N/A	N/A	N/A	
URS149	<p>Move-in and erection Plan 设备进场和安装计划</p> <p>Before system to be delivered. The plan shall be provided and approved. 在系统发货前，设备进场和安装计划应提前递交给业主并得到批准。</p>	I	N/A	N/A	N/A	
URS150	<p>All system has to be accepted by a Factory Acceptance Test (FAT) executed in the supplier's workshop or factory before delivery. 所有系统发货前都应在供应商的生产车间内接受工厂测试 (FAT)</p>	I	N/A	N/A	N/A	
URS151	<p>FAT protocols are to be prepared by the supplier and approved by owner prior to the FAT. The supplier shall execute tests at their premises and owner will indicate which tests will be witnessed by them. FAT 测试文件由供应商准备并在 FAT 前交由业主批准。供应商在工厂进行测试，业主会指出哪些测试需</p>	I	N/A	N/A	N/A	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
	要业主在场。					
URS152	The system will only be released for delivery when all FAT documents are closed out by the owner. 系统只有在业主确认 FAT 阶段可以关闭后再发货。	I	N/A	N/A	N/A	
URS153	Following are some work need to be carried out before FAT: 在主工厂测试 (FAT) 前, 供应商必须完成下列项目: 100% Drawings checks 100% 图纸检查 P&ID compliance check P&ID 图纸合规性检查 GA compliance check 系统布置合规性检查 Pneumatic/vacuum testing of vessels & pipeworks. 气动和真空管路测试 Material of construction verification. 结构材料验证 Surface finish verification. Check on condition of process surfaces including planarity, surface roughness, and drainability etc. 表面处理验证。检查工艺设备表面包括平面, 表面粗糙度和管路排水能力 Slope verification. 坡度确认 All items tagged and labelled as per specification 根据要求提供所有物件的位号和标签 Documentation check; all documents 所有文档检查 All supplier drawings checked and revised to 'for FAT. 检查所有供应商的图纸, 确认后更改为“竣工”版	C	N/A	N/A	will be checked during FAT 将在 FAT 时核查	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS154	During the FAT, a fault log and comment/deviation log must be issued documenting any issues that have been raised/ observed during FAT on a daily basis. FAT 过程中应准备故障记录, 意见/偏差日志以记录每天 FAT 测试时发现的问题。	I	N/A	N/A	N/A	
URS155	The FAT shall include, but not be limited to, the following tests: FAT 测试应包含, 但不限于下列各项: Fabrication & assembly inspection 生产和装配检查 Including mechanical verification and dimension verification to approved drawings, P&ID walk down and components inspection, Insulation and labelling, Welding, Slope verification, Dead-leg absence check in process piping, Drainability, Installation orientation etc. 与批准图纸核对设备及尺寸, P&ID 走向和部件检查, 保温和标识, 焊接, 倾角检查, 工艺管路死角检查, 排水能力, 管路安装方向等。 Riboflavin spray coverage test 核黄素覆盖试验 Components operation performance 部件运行性能测试 Documents, certificates, specifications verification 文件, 证书, 规范确认 Before FAT execution, the pre-requisites of FAT should be satisfied, including but not limited to : The material certificate, welding record, cleaning and passivation record, etc. which belongs to the tank and tube system available on supplier's site. 在 FAT 执行前, FAT 的前提条件必须具备, 包括但不限于: 生产现场测试管罐及仪表部件的材质证书、焊接记录, 清洗钝化相关测试记录等。	C	N/A	N/A	will be checked during FAT 将在 FAT 时核查	
URS156	In this project, module mechanical part and automation part are in charged by different suppliers. The buffer preparation and holding module supplier will take responsibility for the mechanical design, procurement and fabrication as well as the instrumentation in the scope, while the owner pointed dedicated automation	I	N/A	N/A	N/A	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
	system supplier will take responsibility for process control system design. 在此项目中，机械模块和自控部分是由不同供应商负责。缓冲液配制和储液模块供应商负责机械设计，采购和制造以及相关仪表；自控供应商负责控制系统设计。					
URS157	During module design & fabrication phase, the buffer preparation and holding module supplier shall reserve space in module for the control cabinets, power cabinets and HMI etc. which are supplied by automation system supplier. Meantime, the buffer preparation and holding module supplier shall install cable tray and/or cable conduits in the suitable & accessible place within the module frame for further step electrical / pneumatic wiring by automation system supplier which will be defined after contractual condition fixed. 在模块设计&制造阶段，缓冲罐和储液罐供应商应在控制柜，配电柜和 HMI 内留出给自控系统供应商的空间。同时，缓冲罐和储液罐供应商应该安装电缆槽和/或导管。	C	N/A	N/A	will be checked during FAT, SAT 将在 FAT,SAT 时核查	
URS158	Respectively, both supplier shall be completed the design, procurement, installation and testing within their own scope of work. The automation system supplier will send the control cabinet & power cabinet to buffer preparation and holding module system supplier's factory for integration. The tested cabinets will be installed by buffer system supplier at the reserved space in module, at the witness by the automation system supplier. Then the automation supplier shall complete wiring work to the valves, instrumentations or motors within the restriction of cable tray or cable conduit which have already be installed by buffer system supplier. 分开来讲，两家供应商都需完成各自合同范围内的设计、采购、安装、测试工作。自控供应商将控制柜、电机柜发货至设备供应商工厂进行整合。测试柜在自控供应商现场监督下，有设备供应商负责安装在预留位置。自控供应商在设备供应商已经安装好的电缆桥架内完成阀门、仪表、电机的布线工作。	C	N/A	N/A	will be checked during FAT, SAT 将在 FAT,SAT 时核查	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS159	<p>A pre-commissioning with temporary electricity and utility provision shall be done in buffer system supplier's factory before delivery. An Integrated Commissioning Manual including all checklists, test methods and procedures shall be generated by owner's commissioning team by integrating the submission from both suppliers. The commissioning plan shall be reviewed and approved prior to execution. The equipment operation manual shall be provided before commissioning.</p> <p>供应有临时电力和公用工程的预调试在交付前已经在设备供应商工厂内完成。</p> <p>调试执行前调试计划需要得到审核和批准。调试前必须提供设备操作手册。</p> <p>This part shall be integrated by the tank supplier. Because self-control are tube tank suppliers, in M2 project, the owner not to integrate the content of commissioning and validation, also won't bear before M1b project due to hardware, control the joint test that resulted from different vendors, this section in M2 project shall have tube tank testing/verification department to responsible for performing, owner supervision.</p> <p>此部分应由管罐供应商整合。</p> <p>由于自控属于管罐供应商的供应商，在 M2 项目中，业主不会来整合调试和确认的内容，也不会再承担之前 M1b 项目中由于硬件、自控来自不同供应商而导致的联合测试，此部分在 M2 项目中均应有管罐的调试/验证部门来负责执行，业主监督。</p>	I	N/A	N/A	will be checked during FAT, SAT 将在 FAT,SAT 时核查	
URS160	<p>Within one month after the installation. The system SAT shall be executed by supplier and owner will be in attendance.</p> <p>设备安装完成的 1 个月内。开始 SAT 工作时，业主和供应商都将在场。</p>	I	N/A	N/A	N/A	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS161	SAT test protocols are to be prepared by the supplier and approved by owner prior to the SAT. SAT 测试工作文件有供应商准备，业主批准 SAT 文件后才开始 SAT。	I	N/A	N/A	N/A	
URS162	Final and successful SAT report will be provided by supplier. 最终版 SAT 文件和 SAT 报告由供应商提供。	I	N/A	N/A	N/A	
URS163	In general FAT tests shall be repeated. Some tests will be confirmed by spot check. 通常情况 FAT 的测试内容可能会在 SAT 阶段重复，一些测试内容需要现场确认。	C	N/A	N/A	Will be checked during SAT 将在 SAT 时核查	
URS164	Protocols for commissioning and the site acceptance test (SAT) shall be produced by the supplier and approved by owner. 调试和 SAT 文件有供应商提供，并得到业主批准。	I	N/A	N/A	N/A	
URS165	Instruments calibration at site. 现场仪表校验	C	N/A	N/A	Will be calibrated at site, and certificates are required 将在 IQ 前完成仪表现场校验，并提供校验证明	
URS166	IQ/OQ protocols will be prepared by supplier with owner approval and the test execution is by supplier with witness by owner. The supplier shall be responsible for troubleshooting within scope of work. Validation procedure and timeline refer to the VMP provided by owner at a later date. 检查文件由供应商准备并经业主审核批准。 验证程序和时间参考业主稍晚时间提供的验证主计划（VMP）。 IOQ work, including phase water test, shall be completed by the verification team established by the pipe tank supplier, including scheme and execution, which can be assisted by the owner personnel.	I	N/A	N/A	N/A	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
	IOQ 工作包括 phase 水试均应由管罐供应商组建的验证团队负责完成，包括方案及执行，业主人员可以辅助执行。					
URS167	Both buffer preparation and holding module suppliers and the owner appointed automation supplier shall be responsible for executing tests during all stages of commissioning or acceptance activities, with witness as required from the owner team. 缓冲配制和储液模块供应商和业主选择的自控供应商，负责所有阶段的调试和验收活动，业主需要见证。	I	N/A	N/A	N/A	
URS168	Both suppliers shall take responsibility on communication and coordination with each other during design, fabrication, commissioning and qualification period for the right interface information exchange. 双方供应商都有责任在设计，制造，调试和确认阶段沟通和协调正确的信息交流。	I	N/A	N/A	N/A	
URS169	Spare parts and accessories for 2-year operation with detailed list (part no., amount, ISO no. (if applicable), price and delivery time). 备品备件应足够两年使用并附上备品清单（清单包括备件代号，数量，ISO 号(若有)，价格和供货周期）。	C	N/A	N/A	Provide spare parts list in deliverable package. 将在竣工前检查	
URS170	The supplier shall take responsibility for all handling and lifting operations and physical assembly of process equipment components. The supplier shall be responsible for performing checks and tests on the equipment to confirm correct installation. 供应商应监督设备的搬运和起重操作，以及各部件的组装。供应商还有责任进行设备的性能检查和测试，以确定设备安装正确。	I	N/A	N/A	N/A	
URS171	The supplier shall cooperate in planning the installation work. 供应商应协助制定设备安装计划。	I	N/A	N/A	N/A	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS172	The following work will be performed by Others: 下述工作由其他方面完成: General construction works including civil, structural and architectural. 一般的建造工作包括土建, 结构和建筑方面	I	N/A	N/A	N/A	
URS173	In the design of the equipment or system, the position of the weight during the calibration of the weighing module should be considered to facilitate the operation without damaging the tank, components and instruments. 设备或系统设计时应考虑称重模块校验时砝码的放置位置, 方便操作, 且不伤害罐体、部件和仪表。	C	容 器 图 纸 Vessel Draw	WC02-26314 WR02-26315 WR02-26316 WR02-26328	The supporting leg of the weighing equipment is provided with a tensile rod, which makes it easy to place the tray of weighing weights and avoids the meters and other relevant parts 带称重的设备支腿加了拉筋杆, 便于放置称重砝码的托盘, 并避开了仪表及其他相关部件	
URS174	All parts of the equipment must be readily accessible for maintenance, servicing and cleaning. The supplier shall provide lifting points for equipment that must be removed for routine maintenance. The supplier shall advise where lifting aids are required to maintain equipment. 设备的所有部件必须易于维护, 维修和清洗。供应商应为日常维护中需要移动的部件提供起重设备。供应商应告知哪些设备在维护时需要起重设备。	C	Equipment Layout 设备布置图 Piping Layout 管道布置图	PP20-0098-T-009913~009914 PP20-0098-T-009905~009906 PP20-0098-T-009916~009921 PP20-0098-T-009923	All parts of the equipment are readily accessible for maintenance, servicing and cleaning. SMP provide lifting points for equipment that must be removed for routine maintenance. 设备的所有部件必须易于维护, 维修和清洗。供应商应为日常维护中需要移动的部件提供起重设备。	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS175	The supplier's maintenance manual shall include a preventive maintenance plan with times, frequencies and costs for each activity. 供应商提供的维护手册应包含预防性维护计划和周期，频率和每次维护费用。	C	N/A	N/A	Will be checked before handover 将在竣工前核查	
URS176	The supplier shall provide drawings to show access requirements to operate and maintain the equipment. 供应商应提供图例告知操作人员操作和维护设备的要求。	C	N/A	N/A	Will be checked before handover 将在竣工前核查	
URS177	The supplier shall provide details of support available either in or close to site. 供应商应告知在现场附近可提供的详细辅助工作。	I	N/A	N/A	N/A	
URS178	The supplier shall include in proposal the cost of adequate training of operators, supervisors and maintenance staff. 供应商应在其提交的建议书里包含培训业主方操作人员，监管人员和维修人员的费用。	I	N/A	N/A	N/A	
URS179	The supplier shall guarantee that the equipment part of the scope of supply meet design and performance requirements specified, and alter and/or replace, at his own cost, any piece of equipment which fails to meet these requirements (work and factory trained supervision necessary included). 供应商必须保证供货范围内设备部分符合已注明的设计和性能要求。任何部件不能符合上述要求的，由供应商全额负责更换费用（需要有训练过的现场监管人员确认）。	I	N/A	N/A	N/A	
URS180	Warrant all materials and labour included here to be free from defects for a period of 24 months from date of owner. Replace any parts found defective due to manufacture and reinstall new ones at no cost to owner. The warranty period shall be calculated from the date when PQ of equipment or system is completed. 保证与业主商定的时间起 24 个月内，设备的材料和劳务不会出现缺陷。替换任何发现有缺陷的部件，由	I	N/A	N/A	N/A	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
	此造成制造和安装费用与业主无关。 质保期自设备或系统完成 PQ 之日起计算。					
URS181	The type of packing to be used shall be selected with due regard to the chosen shipping media to be used and the environmental conditions to be encountered during shipping, handling and storage. 包装类型的使用应充分考虑到所选择的航运方式，环境因素和在运输，装卸和储存过程中可能遇到的问题。	I	N/A	N/A	N/A	
URS182	Equipment shall be adequately protected during shipping to site. All loose parts shall be adequately boxed crated or bagged. 设备在运输至现场过程中已应充分保护。所有松散的零部件应考虑装盒，装袋或装箱。	I	N/A	N/A	N/A	
URS183	Small items such as bolts, nuts washers, shims, packers and small items of equipment shall be provided in waterproof grit free containers. 如螺栓，螺母垫圈，垫片，设备部件等小件物品应提供无沙防水容器。	I	N/A	N/A	N/A	
URS184	All rotating equipment such as motors, fans etc., in which moving parts could be damaged due to shipping vibration must be secured per the manufacturer's recommendations. 所有可能转动的设备如电机，风机等，其中运动部件会由于运输振动而损坏，应根据生产厂商建议进行妥当保护。	I	N/A	N/A	N/A	
URS185	Pipe, tubing and all openings shall be end capped to prevent the ingress of moisture, dust and any foreign matter that may contaminate the process fluid. Pipe fittings and smaller equipment which contacts the process fluid shall be bagged and tagged to prevent contact with moisture and dust. 管道和管型部件应有封头防止湿气，灰尘等外界异物影响管道质量。可能会接触工艺流体的管件和小部件应袋装保护并有标识，防止其接触湿气和灰尘。	I	N/A	N/A	N/A	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS186	Exposed machined and/or polished surfaces shall be protected with a strippable membrane or coated with a suitable protective compound that shall be easily removable without the use of solvents. 暴露的机械加工和抛光表面可以用易剥落的膜或不用试剂就可除去的化合物涂层保护。	I	N/A	N/A	N/A	
URS187	The method of delivery of equipment associated with this contract is to be organised with owner. 本合同中相关设备的发货方法由业主组织。	I	N/A	N/A	N/A	
URS188	The cost of delivery of the equipment is within the supplier's scope of supply. 设备运输费用在供应商供货范围内。	I	N/A	N/A	N/A	
URS189	The cost of insurance of the equipment delivery is within the supplier's scope of supply. 设备运输途中的保险费用在供应商供货范围内。	I	N/A	N/A	N/A	
URS190	The supplier shall provide complete documents showing and describing in detail the equipment being furnished, including mechanical, construction and equipment detail, and other pertinent data as per design documents, for checking and approval. 供应商应提供完整的文件来详细展示设备生产装配过程，包括机械，结构和设备细节，以及其他相关设计数据，以备检查和批准。	C	N/A	N/A	will be checked during FAT, SAT 将在 FAT,SAT 时核查	
URS191	The supplier documentation shall contain sufficient information to enable owner to proceed with engineering activities. 供应商文件应包含足够的信息保证业主工程活动能正常开展。	C	N/A	N/A	will be checked during FAT, SAT 将在 FAT,SAT 时核查	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS192	Number of copies in English: 英文版文件副本数: Three (3) hard copies as "AS BUILT" documentation plus electronic form (CD). Drawings and documents shall be made in AUTOCAD and Microsoft Office. The Sub-supplier documentation that will be not available in electronic form shall be scanned. 3 分纸质文件作为“竣工”文档，另加电子版（CD）图纸和文件分别以 AUTOCAD 和 Microsoft Office 2007 格式。分包供应商不能提供电子版文档的，应提供扫描版。	C	N/A	N/A	Will be checked before handover 将在竣工前核查	
URS193	Three (3) hard copy of Installation, Operating and Maintenance Manuals, Certification Dossier and Technical Manual in English shall be made available at site before start-up for comments/approval. 3 份纸质文件，安装，运行和维护手册，证书和技术手册，必须在现场开机前审批。	C	N/A	N/A	Will be checked before handover 将在竣工前核查	
URS194	The supplier documents shall be written in both Chinese and English (Chinese mainly) and all units of measurement shall be metric units. 供应商文件以中英文编写（中文为主），所有测量单位为公制单位。	C	N/A	N/A	will be checked during FAT, SAT 将在 FAT,SAT 时核查	
URS195	The supplier shall adopt only his own numbering system of the drawings. 供应商采用自己的编号系统。 The supplier shall adopt only his own numbering system for lines and components. 供应商在提供的管线和组件上采用自己的编号系统。 Equipment number 设备编号 Supplier job number 供应商作业号 Supplier document number and revision. 供应商文件编号和版本号 The documents will not be acceptable unless all required information is incorporated. 文件信息必须完整，才能被接受。	C	N/A	N/A	will be checked during FAT, SAT 将在 FAT,SAT 时核查	

URS 编号 URS No.	URS 描述 URS Description	URS 分类 URS Classification (Q, C, I)	设计文件名称 Design Doc. Name	设计文件编号 Design Doc. No.	设计描述 Design Description	是否可接受 Acceptable Yes 是/No 否
URS196	The supplier shall provide a documented response to any technical and/or commercial issues within one (1) week from receipt of inquiry from owner. 供应商应可以在一周内答复对于来自业主关于技术或商业问题的询问。响应可记录。	I	N/A	N/A	N/A	
URS197	The supplier is not relieved from his responsibility for proper detailing of the design specified on owner contract documents by approval of owner documents. 供应商有责任提供合适的详细设计说明文件直至在业主批准。	I	N/A	N/A	N/A	
URS198	The supplier is not relieved from his responsibility for satisfactory construction, compliance with the design documents and applicable codes, for errors of omissions of any kind in the final product by approval of owner of documents. 供应商有责任提供符合要求的施工，符合文件和应用规范的设计，以及任何需要提交的最终文件中遗漏的部分，直至业主批准。	I	N/A	N/A	N/A	

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可接受的标准 Acceptance Criteria		是 Yes /否 No	
URS 中的要求都有设计文件进行响应，且可接受。 The requirements in URS have been responded to by design documents, and the response is acceptable.			
备注 Comments:			
偏差编号 Deviation No.			
执行人 Executed by		日期 Date	
确认人 Verified by		日期 Date	

10 附件清单 APPENDIX LIST

附件编号 Appendix No.	附件描述 Appendix Description	总页数 Total page

11 执行结果的审批 EXECUTION REVIEW AND APPROVAL

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The Design Qualification activities has been executed completely, and all of DQ checklists have been fulfilled. The results of tests are summarized in the table below:

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Based on the above execution results, the conclusion of the DQ is as follows:

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The execution and the results fulfill the requirements. No any open deviation is present. The system has been authorized to proceed to next qualification step.

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The execution and the results do not completely fulfill the requirements. Open deviations are present, but do not affect the final result of the qualification. The equipment has been authorized to proceed to next qualification step.

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The execution and the results do not fulfill the requirements. Open deviations are present, and do affect the final result of the qualification. The equipment has been not authorized to proceed to next qualification step. Further measures have to take place. Corrections have to be documented separately.

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The signatures below indicate that the DQ has been completed.

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Shanghai Morimatsu Pharmaceutical Equipment Engineering Co., Ltd.

	职务 Function	签名 Signature	日期 Date
审核 Reviewed by	验证经理 Validation Manager		
审核 Reviewed by	设计经理 Design Manager		

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Innovent Biologics (Suzhou) Co., Ltd.

	职务 Function	签名 Signature	日期 Date
审核 Reviewed by			
审核 Reviewed by			
审核 Reviewed by			
审核 Reviewed by			
批准 Approved by			