



sartorius stedim  
biotech

## Sartocheck® 4 | Sartocheck® 4 plus Automatic Detection of Improper Test Setup



Application  
Note

Filter Integrity  
Testing

turning science into solutions

Dr. Jens Meyer  
 Gary Dawson  
 Martin Dahlberg

## Background Information

Filter integrity testing is a crucial requirement in biopharmaceutical production processes. Automatic test units perform indirect tests to verify the pore size to be as specified for a sterilizing-grade membrane filter. Different test methods are established for which a specific test pressure is applied to a wetted membrane. In most cases self-closing quick connectors (e.g. Stäubli couplings) are used for pneumatic connection of the test tubing from the integrity tester to the test housing. In case of an inadvertently disconnected filter housing, the normally-closed coupling will cause the test unit to only pressurize the internal test tubing and not the filter housing. In that situation, the integrity tester will not detect any significant diffusion rate or any relevant nonlinearity in flow, which is needed to detect the bubble point. However, as the only criteria for the test evaluation, a maximum diffusion rate and | or minimum bubble point is defined. Both criteria (max. diffusion rate is not exceeded and bubble point is greater than minimum BP) might be fulfilled when only the tubing is tested. In that situation an integrity test without any filter might be classified as "test passed," which indeed is a false positive test result.

An experienced operator may immediately recognize the "flatline" of the BP or diffusion curve. Nevertheless, there is a remaining risk for the daily routine work. This application note describes a new feature of the Sartocheck® 4 and Sartocheck® 4 plus which increases the detection capability of the integrity testing procedure by identifying an improper test setup.

## How to Avoid False Positive Test Results?

In contrast to previous versions, the latest software version of Sartocheck® 4 (version 2.03) and Sartocheck® 4 plus has added three independent parameters that provide the enhanced detection of an incorrect test set up. The proper use of these parameters will reliably detect situations in which the test set up is not pneumatically connected to the integrity tester, or if a blind end cap or pinched line is preventing normal air flow.





Figure 1:  
 Screenshot of Sartocheck® 4 that shows how to activate the additional parameters: a non-standard test has to be chosen.

To preserve the well established foundation of the test pass | fail status, the criteria used to define the outcome of the tests have not been changed. However, when an improper test setup is detected, an additional warning message will subsequently be displayed and printed in red with the



following text:

Figure 2:  
 Screenshot of Sartocheck® 4 showing the warning message that is displayed if an incorrect test setup is detected.

The following parameters to detect an improper test set up can be entered by the user:

**Min. Diffusion**

(for the Diffusion, WIT and WFT tests)  
This parameter defines a realistic minimum diffusion rate that would be expected for a filter that is properly connected to the Sartocheck® (there must be some amount of expected diffusion if a filter is connected).

**Min. Net Volume** (for all tests that include volume determination step)

This parameter defines the minimum upstream volume of the test setup (tubing plus filter capsule or filter housing). A test set up connected to the Sartocheck® via the test tubing will have a specific net upstream volume. A disconnected filter (or filter housing) will result in an inordinately low upstream volume that only reflects the internal tubing of the Sartocheck unit itself.

**Min. Flow BP Test**

(for the Bubble Point test only)  
This parameter defines the diffusive flow rate just prior to reaching BP<sub>max</sub>. A disconnected filter (or filter housing) will result in an inordinately low flow rate.

The concept behind these three parameters is that as soon as any filter housing is connected to the Sartocheck, and a test is started, an increased net volume and a significant diffusion rate should be measured. To cover all possible conditions, all three parameters are required.




Figure 3:  
Screenshot of the Sartocheck® 4 (version 2.03) which shows where the additional parameters are defined.

**Adjusting the parameters**

These three parameters have been implemented to detect a situation in which a filter (housing or capsule) has not been properly connected to the Sartocheck®. As such, these parameters should be programmed to distinguish between the situations of a “filter connected” and “no filter connected.” Because the respective flow rates and net volumes are dependent upon on the individual test setup, no predetermined values can be recommended. Rather, the general recommendations for selecting the optimal values are as follows:

**Minimum Diffusion:**  
Acquire the “blind diffusion” (in mL/min.) by running a diffusion test with no test setup connected. This measured diffusion can be used as Minimum Diffusion value. This value should always be very low (something around 3 or 4 mL/min.) and should never be near the true diffusion of a passing filter cartridge to avoid false error messages.

**Minimum Flow BPT:**  
The parameter min diffusion BP test is the diffusion rate that is measured just before reaching the Max. BP pressure. As such, a good rule for defining this value is using 1/2 of the Diff<sub>max</sub> value.

**Minimum Net Volume:**  
This value should be programmed within 80–90% of the actual net upstream volume of the housing with the filter installed, or the upstream volume of the filter capsule (if no housing is used). For example: If the SC4 measures a net upstream volume of 1000 mL for a given test setup, then a Minimum Net Volume of 850 mL should be programmed. The only requirement is to define a value that lies between the respective volumes for a decoupled test set up, and the net volume of a properly connected test set up, respectively. This is typically a very broad range, making the Minimum Net Volume very easy to define somewhere in between the two extremes.

The Sartocheck® 4 allows the choice of either activating or deactivating the new safety parameters. However, in case of a decoupled test setup, this condition will only be reliably identified when the parameters are activated:




Figure 4:  
Printouts of Sartocheck® 4 with deactivated (left) and activated (right) safety parameters. A diffusion test was performed with a disconnected test tubing. This condition was only identified when the new parameters were activated, resulting in a warning message displayed on the screen and printed in red.

**Conclusion:**

The new detection features of Sartocheck® 4 (≥ software version 2.03) and Sartocheck® 4 plus significantly increase the robustness of filter integrity testing. When those additional three parameters (Min. Diffusion, Min Flow BPT, and Min Net Volume) are defined correctly, the Sartocheck® reliably eliminates the possibility of false positive test results that might otherwise occur in the daily work routine due to a disconnected or improper test setup.

# Sales and Service Contacts

For further contacts, visit [www.sartorius-stedim.com](http://www.sartorius-stedim.com)

## Europe

### Germany

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen

Phone +49.551.308.0  
Fax +49.551.308.3289

[www.sartorius-stedim.com](http://www.sartorius-stedim.com)

Sartorius Stedim Systems GmbH  
Schwarzenberger Weg 73-79  
34212 Melsungen

Phone +49.5661.71.3400  
Fax +49.5661.71.3702

[www.sartorius-stedim.com](http://www.sartorius-stedim.com)

### France

Sartorius Stedim Biotech S.A.  
ZI Les Paluds  
Avenue de Jouques – BP 1051  
13781 Aubagne Cedex

Phone +33.442.845600  
Fax +33.442.845619

Sartorius Stedim France SAS  
ZI Les Paluds  
Avenue de Jouques – CS 71058  
13781 Aubagne Cedex

Phone +33.442.845600  
Fax +33.442.846545

### Austria

Sartorius Stedim Austria GmbH  
Franzosengraben 12  
A-1030 Vienna

Phone +43.1.7965763.18  
Fax +43.1.796576344

### Belgium

Sartorius Stedim Belgium N.V.  
Leuvensesteenweg, 248/B  
1800 Vilvoorde

Phone +32.2.756.06.80  
Fax +32.2.756.06.81

### Denmark

Sartorius Stedim Nordic A/S  
Hoerskaetten 6D, 1.  
DK-2630 Taastrup

Phone +45.7023.4400  
Fax +45.4630.4030

### Italy

Sartorius Stedim Italy S.p.A.  
Via dell'Antella, 76/A  
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41  
Fax +39.055.63.40.526

### Netherlands

Sartorius Stedim Netherlands B.V.  
Edisonbaan 24  
3439 MN Nieuwegein

Phone +31.30.6025080  
Fax +31.30.6025099

### Spain

Sartorius Stedim Spain SA  
C/Isabel Colbrand 10,  
Oficina 70  
Polígono Industrial de Fuencarral  
28050 Madrid

Phone +34.90.2110935  
Fax +34.91.3589623

### Switzerland

Sartorius Stedim Switzerland GmbH  
Lerzenstrasse 21  
8953 Dietikon

Phone +41.44.741.05.00  
Fax +41.44.741.05.09

### U.K.

Sartorius Stedim UK Limited  
Longmead Business Park  
Blenheim Road, Epsom  
Surrey KT19 9 QQ

Phone +44.1372.737159  
Fax +44.1372.726171

## America

### USA

Sartorius Stedim North America Inc.  
5 Orville Drive  
Bohemia, NY 11716

Toll-Free +1.800.368.7178  
Fax +1.631.254.4253

Sartorius Stedim SUS Inc.  
1910 Mark Court  
Concord, CA 94520

Phone +1.925.689.6650  
Toll Free +1.800.914.6644  
Fax +1.925.689.6988

Sartorius Stedim Systems Inc.  
201 South Ingram Mill Road  
Springfield, MO 65802

Phone +1.417.873.9636  
Fax +1.417.873.9275

### Argentina

Sartorius Argentina S.A.  
Int. A. Avalos 4251  
B1605ECS Munro  
Buenos Aires

Phone +54.11.4721.0505  
Fax +54.11.4762.2333

### Brazil

Sartorius do Brasil Ltda  
Av. Dom Pedro I, 241  
Bairro Vila Pires  
Santo André  
São Paulo  
Cep 09110-001

Phone +55.11.4451.6226  
Fax +55.11.4451.4369

### Mexico

Sartorius de México S.A. de C.V.  
Circuito Circunvalación Poniente No. 149  
Ciudad Satélite  
53100 Naucalpan, Estado de México

Phone +52.5555.62.1102  
Fax +52.5555.62.2942

## Asia | Pacific

### China

Sartorius Stedim Beijing  
Representative Office  
No. 33, Yu'an Road,  
Airport Industrial Zone B, Shunyi District  
Beijing 101300

Phone +86.10.80426516  
Fax +86.10.80426580

Sartorius Stedim Shanghai  
Representative Office  
Room 618, Tower 1, German Centre,  
Shanghai, PRC., 201203

Phone +86.21.28986393  
Fax +86.21.28986392.11

Sartorius Stedim Guangzhou Office  
Room 704, Broadway Plaza,  
No. 233-234 Dong Feng West Road  
Guangzhou 510180

Phone +86.20.8351.7921  
Fax +86.20.8351.7931

### India

Sartorius Stedim India Pvt. Ltd.  
10, 6th Main, 3rd Phase Peenya  
KIADB Industrial Area  
Bangalore – 560 058

Phone +91.80.2839.1963|0461  
Fax +91.80.2839.8262

### Japan

Sartorius Stedim Japan K.K.  
KY Building, 8-11  
Kita Shinagawa 1-chome  
Shinagawa-ku  
Tokyo 140-0001

Phone +81.3.3740.5407  
Fax +81.3.3740.5406

### Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.  
Lot L3-E-3B, Enterprise 4  
Technology Park Malaysia  
Bukit Jalil  
57000 Kuala Lumpur

Phone +60.3.8996.0622  
Fax +60.3.8996.0755

### Singapore

Sartorius Stedim Singapore Pte. Ltd.  
10, Science Park Road, The Alpha  
#02-25, Singapore Science Park 2  
Singapore 117684

Phone +65.6872.3966  
Fax +65.6778.2494

### Australia

Sartorius Stedim Australia Pty. Ltd.  
Unit 5, 7-11 Rodeo Drive  
Dandenong South Vic 3175

Phone +61.3.8762.1800  
Fax +61.3.8762.1828