



Coloured P&ID for TFF System Rev1.0

Based on P&ID Drawing 05.000976-PID Rev1.1



CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

Revision History

Revision	Description	Date
1.0	First issue	22 Oct 2020

Phase list

Phase 1: Drain Vessel
Phase 2: Fill vessel with media
Phase 3: Re-circulation filtration loop
Phase 4: Re-circulation vessel spray ball
Phase 5: CIP Inlets
Phase 6: Blow down pipe work
Phase 7: Flush membrane cassettes
Phase 8: Temperature control
Phase 9: NWP-Test
Phase 10: Concentration
Phase 11: Diafiltration with tank level control.
Phase 12: Product Recovery
Phase 13: Retentate re-circulation
Phase 14: System pressure hold test
Phase 15: Integrity test filter cassettes
Phase 16: Concentration with tank level control.
Phase 17: Tank SIP




Phase 1: Drain vessel



CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1Na



Phase	Drain vessel	
Step	Description	Transition
10	DRAIN VESSEL	Empty signal at LS-101

CONNECTIONS




Phase 2: Fill vessel with media



CONFIDENTIAL


This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1N



CONNECTIONS

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1N



CONNECTIONS



Phase 3: Re-circulation filtration loop




CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1N NaOH

B002 = 0 (2 cassette holders)


B002 = 1 (4 cassette holders)



Phase	Re-circulation filtration loop	
Step	Description	Transition
10	Prime Loop	*Timer Z105 expired

CONNECTIONS

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1Na






Phase 4: Re-circulation vessel spray ball



CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Separate Cooling water use
- Add XV403
- Exchange WFI and 0.1N






Phase 5: CIP Inlet



CONFIDENTIAL


This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Separate Cooling Water use
- Add XV403
- Exchange WFI and 0.1Na



CONNECTIONS


- Separate Cooling Water use
- Add XV403
- Exchange WFI and 0.1N



Phase	CIP inlet	
Step	Description	Transition
30	Inlet SIP	*Timer Z102 expired

CONNECTIONS

- Separate Cooling Water use
- Add XV403
- Exchange WFI and 0.1N



CONNECTIONS




Phase 6: Blow-down pipe work



CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.


- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1Na



Phase	Blow-down pipe work	
Step	Description	Transition
10	Pressurize vessel	Pressure at PT501 ≥ Set point P000


CONNECTIONS

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1Na




CONNECTIONS

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1Na



CONNECTIONS

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1Na



Phase	Blow-down pipe work	
Step	Description	Transition
40	Vent system	•Pressure PT501 ≤ Z002

CONNECTIONS



Phase 7: Flush membrane cassettes




CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1Na

B000 = 0 (2 cassette holders)

B000 = 1 (4 cassette holders)



Phase	Flush membrane cassettes	
Step	Description	Transition
10	Flush Retentate	WIT501 actual ≤ Starting volume (V000) – V001

CONNECTIONS

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1N NaOH

B000 = 0 (2 cassette holders)

B000 = 1 (4 cassette holders)

B001 = 0

B001 = 1

Phase	Flush membrane cassettes	
Step	Description	Transition
20 25	Flush Permeate	WIT501 actual \leq Starting volume (V000) - V001 - V002

CONNECTIONS




Phase 8: Temperature control




CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- operator cooling water use
- Add XV403
- Exchange WFI and 0.1Na




- operator cooling water use
- Add XV403
- Exchange WFI and 0.1Na



CONNECTIONS

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1N



CONNECTIONS



Phase 9: NWP test




CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Separate Cooling Water
- Add XV403
- Exchange WFI and 0.1N

B000 = 0 (2 cassette holders)

B000 = 1 (4 cassette holders)



CONNECTIONS

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1Na

B000 = 0 (2 cassette holders)

B000 = 1 (4 cassette holders)

Phase	NWP Test	
Step	Description	Transition
10	Recirculation / Stabilization	TMP-1 set point -controller deviation low \leq TMP-1 at PCV201 \leq TMP set point + controller deviation high
20	Measurement step / clean water flux calculation (close all valves afterwards)	<ul style="list-style-type: none"> • Timer OP Z001 expired \rightarrow go to step 30. • TMP controller is out of high dev. / low dev. Range \rightarrow go to step 10.

CONNECTIONS



Phase 10: Concentration



CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Separate Cooling Water use
- Add XV403
- Exchange WFI and 0.1Na

B004 = 0 (2 cassette holders)

B004 = 1 (4 cassette holders)

Only if B000=1

Phase	Concentration	
Step	Description	Transition
10	Start Recirculation / Fill Loop	Timer OP Z100 expired If B003 = 1 => go to step 20 If B003 = 0 => go to step 30

• Separate cooling water use
 • Add XV403
 • Exchange WFI and 0.1Na

B004 = 0 (2 cassette holders)

B004 = 1 (4 cassette holders)

B001 = 1/2 => C

B001=3 => F

C/F

Only if B000=1

Phase	Concentration	
Step	Description	Transition
20	Adjust process conditions	<ul style="list-style-type: none"> • Only executed if B003 = 1 • TMP-1 set point –controller deviation low TMP-1 at PCV201 TMP-1 set point + controller deviation high • AND • Flow set point –controller deviation low RP101 at FIT 201 Flow set point + controller deviation high • AND • FIT301 set point –controller deviation low FIT301 at PCV301 FIT set point + controller deviation high (only B009 =2) • AND • TMP-2 set point –controller deviation low TMP-2 at PCV301 TMP-2 set point + controller deviation high (only B009 =1) • => go to step 30 if B005 = 1 and actual value 他04 ≥ OP Q000 => go to step 40

CONNECTIONS

Separate Cooling Water use
 Add XV403
 Exchange WFI and 0.1Na

B004 = 0 (2 cassette holders)

B004 = 1 (4 cassette holders)

B001 = 1/2 => C

B001=3 => F

C/F

Only if B000=1

Phase	Concentration	
Step	Description	Transition
30	Concentration	<p>•• If B008 = 0: $(V103 / B007) \geq ((\text{Actual value WIT501}) + V100 + (B006 * V101))$ => go to step 50.</p> <p>•• If B008 = 1: $V110 \geq ((\text{Actual value WIT501}) + V100 + (B006 * V101))$ => go to step 50.</p> <p>•• Only if B005 = 1 and :actual value AT304 > OP Q000 => go to step 40</p> <p>•• Only if B001=1 and operator confirmation=> go to step 35</p>

CONNECTIONS

- Separate Cooling Water use
- Add XV403
- Exchange WFI and 0.1N

B004 = 0 (2 cassette holders)

B004 = 1 (4 cassette holders)

Only if B000=1

Phase	Concentration	
Step	Description	Transition
35	Concentration with miniloop	<p>•• If B008 = 0: $(V103 / B007) \geq ((\text{Actual value WIT501}) + V100 + (B006 * V101))$ \Rightarrow go to step 50.</p> <p>•• If B008 = 1: $V110 \geq ((\text{Actual value WIT501}) + V100 + (B006 * V101))$ \Rightarrow go to step 50.</p> <p>•• Only if B005 = 1 and :actual value AT304 > OP Q000 \Rightarrow go to step 40</p>

CONNECTIONS

• Separate cooling water use
 • Add XV403
 • Exchange WFI and 0.1Na

B004 = 0 (2 cassette holders)

B004 = 1 (4 cassette holders)

Only if B000=1

Phase	Concentration	
Step	Description	Transition
30	UV Alarm (Fall back recirculation)	Operator selection: Selection 1: •Continue => go to step 30 (ignore UV alarm) Selection 2: •End phase => go to step 50 (End) and set variable B101 = 1 (Variable can be used as a criterium in a jump phase in the recipe).

CONNECTIONS



Phase 11: Diafiltration with tank level control



CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Separate Cooling Water use
- Add XV403
- Exchange WFI and 0.1Na

B005 = 0 (2 cassette holders)

B005 = 1 (4 cassette holders)

Only if B006=1

Phase	Diafiltration with tank level control	
Step	Description	Transition
10	Start Recirculation / Fill Loop	Timer OP Z100 expired If B008 = 1 => go to step 20 If B008 = 0 => go to step 30

CONNECTIONS

• Separate cooling water use
 • Add XV403
 • Exchange WFI and 0.1Na

B005 = 0 (2 cassette holders)

B005 = 1 (4 cassette holders)

B001 = 1/2 => C
 B001=3 => F

C/F

Only if B006=1

Phase	Diafiltration with tank level control	
Step	Description	Transition
20	Adjust process conditions	<ul style="list-style-type: none"> • Only executed if B008 = 1 • TMP-1 set point –controller deviation low \leq TMP-1 at PCV201 \leq TMP-1 set point + controller deviation high • AND • Flow set point –controller deviation low \leq RP101 at FIT 201 \leq Flow set point + controller deviation high • AND • FIT301 set point –controller deviation low \leq FIT301 at PCV301 \leq FIT set point + controller deviation high (only B009 = 2) • AND • TMP-2 set point –controller deviation low \leq TMP-2 at PCV301 \leq TMP-2 set point + controller deviation high (only B009 = 1) • if actual value AT304 > OP Q000 and B007=1 => go to step 90 • else => go to step 30

CONNECTIONS

• Separate cooling water use
 • Add XV403
 • Exchange WFI and 0.1N

B005 = 0 (2 cassette holders)

B005 = 1 (4 cassette holders)

B001 = 1/2 => C

B001=3 => F

C/F

Only if B006=1

Phase	Diafiltration with tank level control	
Step	Description	Transition
30	Diafiltration	<ul style="list-style-type: none"> totalised permeate flow at FITQ301 \geq V001 (only if B003=0) => go to step 70 (phase end) totalised permeate flow FIT301 \geq (V103+actual value WIT501+(V104*B011))*B010 (only if B003 = 1) => go to step 70 (phase end) WIT501 \leq Level control setpoint - Level controller deviation low => go to step 50 WIT501 \geq Level control setpoint + Level Controller deviation high => go to step 40 Only if B007 = 1 and AT304 > OP Q000 => go to step 60

CONNECTIONS

• Separate cooling water use
 • Add XV403
 • Exchange WFI and 0.1Na

B005 = 0 (2 cassette holders)

B005 = 1 (4 cassette holders)

B001 = 1/2 => C

B001=3 => F

C/F

Only if B006=1

Phase	Diafiltration with tank level control	
Step	Description	Transition
40	Concentration without media transfer due to high level warning	Vessel Volume $WIT501 \leq \text{Set point } V002 / V102$ => go to step 30 totalised permeate flow at $FITQ301 \geq V001$ (only if B003=0) => go to step 70 (phase end) totalised permeate flow $FIT301 \geq (V103 + \text{actual value } WIT501 + (V104 * B011)) * B010$ (only if B003 = 1) => go to step 70 (phase end) only if B007 = 1 and $AT304 \geq OP Q000$ => go to step 60

CONNECTIONS

• Separate cooling water use
 • Add XV403
 • Exchange WFI and 0.1N

B005 = 0 (2 cassette holders)

B005 = 1 (4 cassette holders)

B001 = 1/2 => C
 B001=3 => F

C/F


C

Only if B006=1

Phase	Diafiltration with tank level control	
Step	Description	Transition
50	Media transfer without concentration due to low level warning	••Vessel Volume WIT501 ≥ Set point V002 => got to step 30

F

F



CONNECTIONS

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1Na

B005 = 0 (2 cassette holders)

B005 = 1 (4 cassette holders)

Only if B006=1

Phase		
Diafiltration with tank level control		
Step	Description	Transition
30	UV Alarm (Fall back recirculation)	<p>Operator selection:</p> <p>Selection 1: Continue => go to step 20 (only B008=1) (ignore UV alarm) (Supervisor only)</p> <p>Selection 2: Continue => go to step 30 (only B008=0) (ignore UV alarm) (Supervisor only)</p> <p>Selection 3: End phase => go to step 70 (End) but continue recipe (Supervisor only)</p>

CONNECTIONS




Phase 12: Product Recovery



CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.


- Separate Cooling water use
- Add XV403
- Exchange WFI and 0.1Na



Phase	Product Recovery	
Step	Description	Transition
10	Temperature control	• Actual value at TT501 \leq T000

CONNECTIONS

- Upgrade Cooling Water
- Add XV403
- Exchange WFI and 0.1Na



Phase	Product Recovery	
Step	Description	Transition
20	Recovery from vessel	•Empty signal at LS101 for Time Z000

CONNECTIONS

- Separate Cooling Water use
- Add XV403
- Exchange WFI and 0.1N

B001 = 0 (2 cassette holders)

B001 = 1 (4 cassette holders)

Only if B000=1

Phase	Product Recovery	
Step	Description	Transition
30	Recovery from filter cassettes	<ul style="list-style-type: none"> •Timer Z101 expired If B003=0 then go to step 50 If B003=1 then go to step 40



Phase 13: Retentate re-circulation




CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1N

B000 = 0 (2 cassette holders)

B000 = 1 (4 cassette holders)



Phase	Retentate re-circulation	
Step	Description	Transition
10	Buffer Recirculation	•Timer Z000 expired

CONNECTIONS



Phase 14: System pressure hold test



CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Separate cooling water system
- Add XV403
- Exchange WFI and 0.1N NaOH

B000 = 0 (2 cassette holders)

B000 = 1 (4 cassette holders)

System pressure hold test		
Phase	System pressure hold test	
Step	Description	Transition
10	Pressurize system	• Actual Pressure at PT501 \geq P000-P002
20	Stabilize system	Timer Z100 expired

CONNECTIONS

• Upgrade cooling water system
 • Add XV403
 • Exchange WFI and 0.1Na

B000 = 0 (2 cassette holders)

B000 = 1 (4 cassette holders)

Phase	System pressure hold test	
Step	Description	Transition
30	Measurement Step	<ul style="list-style-type: none"> Timer Z000 expired → go to step 40 If pressure at PT501 ≤ P003 – P001 → go to step 50
40	Pressure hold test passed	<ul style="list-style-type: none"> Timer Z101 expired → go to step 60
50	Pressure hold test failed	Operator selection: 1: Repeat. → Go to step 10 2: Continue. → Go to step 60

CONNECTIONS

- Separate Cooling Water use
- Add XV403
- Exchange WFI and 0.1N

B000 = 0 (2 cassette holders)

B000 = 1 (4 cassette holders)

Phase	System pressure hold test	
Step	Description	Transition
60	Vent system	PT501 ≤ GP000

CONNECTIONS




Phase 15: Integrity test filter cassettes



CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.


- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1N




Phase	Integrity test filter cassettes	
Step	Description	Transition
10 If B001 = 2	Pressurize system	Pressure at PT201 \geq OP P000

CONNECTIONS


- Repeat cooling water test
- Add XV403
- Exchange WFI and 0.1N



- Repeat cooling water test
- Add XV403
- Exchange WFI and 0.1N




- Upgrade cooling water system
- Add XV403
- Exchange WFI and 0.1N NaOH




CONNECTIONS

• Upgrade cooling water system
 • Add XV403
 • Exchange WFI and 0.1N




CONNECTIONS

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1N




CONNECTIONS

- Separate cooling water use
- Add XV403
- Exchange WFI and 0.1N



CONNECTIONS

- Upgrade cooling water system
- Add XV403
- Exchange WFI and 0.1N NaOH



CONNECTIONS




Phase 16: Concentration with tank level control



CONFIDENTIAL

This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Upgrade cooling water system
- Add XV403
- Exchange WFI and 0.1Na



Phase	Concentration with tank level control	
Step	Description	Transition
10	Enter total product volume	Operator entry

CONNECTIONS

- Separate Cooling Water use
- Add XV403
- Exchange WFI and 0.1N

B002 = 0 (2 cassette holders)

B002 = 1 (4 cassette holders)

Only if B004=1

Phase	Concentration with tank level control	
Step	Description	Transition
20	Start Recirculation / Fill Loop	<ul style="list-style-type: none"> •Timer Z100 expired If B008 = 0 => go to step 40 If B008 = 1 => go to step 30

CONNECTIONS

• Separate Cooling Water use
 • Add XV403
 • Exchange WFI and 0.1Na

B002 = 0 (2 cassette holders)

B002 = 1 (4 cassette holders)

B001 = 1/2 => C
 B001=3 => F

C/F

Only if B004=1

Phase	Concentration with tank level control	
Step	Description	Transition
30	Adjust process conditions	<ul style="list-style-type: none"> • Only executed if B008 = 1 • TMP-1 set point –controller deviation low \leq TMP-1 at PCV201 \leq TMP-1 set point + controller deviation high • AND • Flow set point –controller deviation low \leq RP101 at FIT 201 \leq Flow set point + controller deviation high • AND • FIT301 set point –controller deviation low \leq FIT301 at PCV301 \leq FIT set point + controller deviation high (only B009 = 2) • AND • TMP-2 set point –controller deviation low \leq TMP-2 at PCV301 \leq TMP-2 set point + controller deviation high (only B009 = 1) <p>=> go to step 40 if actual value QIR-30 \geq OP Q000 and B005=1 => go to step 90</p>

CONNECTIONS

• Separate Cooling Water use
 • Add XV403
 • Exchange WFI and 0.1N

B002 = 0 (2 cassette holders)

B002 = 1 (4 cassette holders)

B001 = 1/2 => C

B001=3 => F

C/F

Only if B004=1

Phase	Concentration with tank level control	
Step	Description	Transition
40	Concentration with product transfer	<ul style="list-style-type: none"> Level Controller deviation low => go to step 60 Level Controller deviation high => go to step 50 If B005 = 1 and AT304 ≥ OP Q000 => go to step 90

CONNECTIONS

• Separate cooling water use
 • Add XV403
 • Exchange WFI and 0.1N

B002 = 0 (2 cassette holders)

B002 = 1 (4 cassette holders)

B001 = 1/2 => C

B001=3 => F

C/F

Only if B004=1

Phase	Concentration with tank level control	
Step	Description	Transition
50	Concentration without media transfer due to high level warning	<ul style="list-style-type: none"> Vessel Volume WIT501 ≤ Set point V000 => got to step 40 If B005 = 1 and AT304 ≥ OP Q000 => go to step 90

CONNECTIONS

• Separate cooling water use
 • Add XV403
 • Exchange WFI and 0.1Na

B002 = 0 (2 cassette holders)

B002 = 1 (4 cassette holders)

B001 = 1/2 => C

B001=3 => F

C/F

Only if B004=1

Phase	Concentration with tank level control	
Step	Description	Transition
60	Media transfer without concentration due to low level warning	<ul style="list-style-type: none"> Level Controller deviation low => go to step 60 Level Controller deviation high => go to step 50 If B005 = 1 and AT304 ≥ OP Q000 => go to step 90

CONNECTIONS

- Separate Cooling Water use
- Add XV403
- Exchange WFI and 0.1N

B002 = 0 (2 cassette holders)

B002 = 1 (4 cassette holders)

Only if B004=1

Concentration with tank level control		
Phase	Concentration with tank level control	
Step	Description	Transition
70	Operator confirmation	Operator selection: Selection 1: End => go to End (step 100) Selection 2: Continue => Go to step 30 Selection 3: Final Concentration => Go to step 80 and set variable V109 to (actual value of FITQ301+Actual value WIT501+V100+(B006*V101))/B007)

CONNECTIONS

• Separate Cooling Water use
 • Add XV403
 • Exchange WFI and 0.1Na

B002 = 0 (2 cassette holders)

B002 = 1 (4 cassette holders)

C/F

B001 = 1/2 => C
 B001=3 => F

Only if B004=1

Phase		
Concentration with tank level control		
Step	Description	Transition
80	Final Concentration without media transfer	<ul style="list-style-type: none"> V109 $\geq ((\text{Actual value WIT501}) + V100 + (B006 * V101)) \Rightarrow$ go to step 100. Only if B005 = 1 and actual value AT304 \geq OP Q000 \Rightarrow go to step 90. <p>Only if B001=1 and operator confirmation \Rightarrow go to step 85</p>

CONNECTIONS

• Separate cooling water use
 • Add XV403
 • Exchange WFI and 0.1Na

B002 = 0 (2 cassette holders)

B002 = 1 (4 cassette holders)

Only if B004=1

Phase		
Concentration with tank level control		
Step	Description	Transition
85 (only B001=1)	Concentration with miniloop	<ul style="list-style-type: none"> $V109 \geq ((\text{Actual value WIT501}) + V100 + (B006 * V101)) \Rightarrow$ go to step 100. Only if B005 = 1 and actual value AT304 \geq OP Q000 \Rightarrow go to step 90. <p>Only if B001=1 and operator confirmation \Rightarrow go to step 85</p>

CONNECTIONS

- Separate Cooling Water use
- Add XV403
- Exchange WFI and 0.1N

B002 = 0 (2 cassette holders)

B002 = 1 (4 cassette holders)

Only if B004=1

Phase		
Concentration with tank level control		
Step	Description	Transition
90	UV Alarm (Fall back recirculation)	<ul style="list-style-type: none"> •Operator selection Selection 1: Continue => go to step 30 (ignore UV alarm) (supervisor only) Selection 2: End phase => go to step 100 (End) but continue recipe (supervisor only)

CONNECTIONS




Phase 17: Tank SIP



CONFIDENTIAL


This presentation is the confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2019 Pall Corporation.

- Upgrade Cooling Water
- Add XV403
- Exchange WFI and 0.1N




CONNECTIONS

- Upgrade Cooling Water
- Add XV403
- Exchange WFI and 0.1N




CONNECTIONS

- Upgrade Cooling water use
- Add XV403
- Exchange WFI and 0.1Na




CONNECTIONS

- Upgrade cooling water use
- Add XV403
- Exchange WFI and 0.1N




Phase	Tank SIP	
Step	Description	Transition
50	SIP cooling	TT101/TT102/TT501<P002

- Upgrade Cooling Water
- Add XV403
- Exchange WFI and 0.1N



CONNECTIONS

- Upgrade Cooling Water
- Add XV403
- Exchange WFI and 0.1N



Phase	Tank SIP	
Step	Description	Transition
70	Final Vent	timer Z104 expired